

OPTOPRIM

LASER

engineering overview

MUSE meeting INFN

12 May 2017

Salvatore Salerno



OPTOPRIM srl

Optoprim is a distributor for optoelettronic components, and laser sub systems. We have offices in France, German and Italy.

Optoprim S.r.l. has been founded in 1999

9 people in Monza (main office)

3 people in Roma

4 senior Sales Engineer; 2 junior Sales Engineer



3 Tecnichian for installation, service and process development

We have an internal lab with laser and equipment to prepare samples and for process validation in accordance with customers requests.



Market and products

Laser Processing

Laser: fiber, DPSS, ps, fs, direct diodes Process head for: cutting, welding and cladding Laser machine for high end custom applications optics, F-theta, galvanometer

Scientific/R&D

Instrumentation (laser and light sources analysis and characterization), Coherent and wide spectrum light sources Non-destructive superficial and volumetric analysis of materials Laser protection systems

Medical:

Diodes: free space ed fiber coupled, optics, protection goggles, power meter

Military, defence and aerospace:

Laser, optics, Instrumentation



Product range:





LASER

Summary



Concepts of laser theory and laser light features.

Main laser types

Principal laser applications



LASER

Light Amplification by Stimulated Emission of Radiation

1917: Einstein states the stimulated emissions theory to explain the Plank black body law

1954: Gordon designs the MASER (Microwave Amplification

1960: Basov, Krokhin e Popov develope the laser theory

1960: Maiman manufacture the first laser (ruby based)







If energy is supplied to the atom (absorption), an electron can "jump" to an higher energy state. The atom gets excited, and has more energy than before. On the opposite, if the electron decays to a lower energy level, the atom decreases its energy and releases a photon.





Spontaneous emission

This process happens randomly: photons (of different energy) are emitted in all directions and the associated waves are out of phase.

For example, in a light bulb the current warms the filament, bringing its atoms into a variety of excited state.

Subsequently, they decay by emitting light.





Stimulated emission



The stimulated emission consists in the decaying of the excited level by a photon of the same energy (theorized by Einstein).

Thus, two quantizes of identical em radiations are obtained, that is, coherent (with the same phase), with the same wavelength and in the same direction.



energy radiation equal to the excitation level, the electrons will immediately fall into the lower orbit, emitting photons of the same energy, phase and direction of the incident one.







Conditions required

To generate a laser from the stimulated emission process, it is required to have an electron accumulation in the excited state of the system. To achieve this status, the following three conditions must be met:

1) Inversion of population;

2) The excited state must be a metastable state;

3) The emitted photon must be able to stimulate other photons of the system and must therefore be confined.





1. Inversion of population

Atoms are usually in the fundamental state, so an incident radiation has a greater chance of being absorbed rather than producing stimulated emission.

It is therefore necessary to carry several atoms in a higher energy state, leaving almost empty the lowest state .

This condition is called POPULATION INVERSION.

An incident photon could cause an "avalanche" of stimulated electrons, all perfectly in phase. Such stimulating wave would continue to increase in intensity while propagating through the active medium, as long as the inversion of population can be maintained.







That is, when the number of atoms excited in the active medium is greater than the number of atoms in the fundamental state.

In a two-level laser design you can get at the most: N2 = N1

This is therefore not a beneficial device.

Lasers structures are then designed as 3 or 4 level systems.





2. metastable excited state

If the excited state is metastable (that is, it has a long average life compared to the short time of excited states $\sim 10^{-8}$ s), then the stimulated emission is more likely to occur than the spontaneous emission.







3. Resonant cavity



One of the mirrors has a R \sim 100%, the other has a R <100% to allow the light to go out as laser beam. The rest remains trapped and continues the stimulated emission process.

The two mirrors, which confine the photons into the cavity, form a resonant cavity: they allow to propagate only wavelengths which are an exact multiple of the cavity length.

In addition, only photons directed along the cavity axis are propagated and amplified.











There are several technologies used to generate pulsed laser, exploiting which today is possible to generate ultra-short pulses (the record is 67 as 10^{-18} s) or amplified systems that emit peak power of PW (10^{15} W)



Laser light properties: Unidirectional

A laser beam propagates with a minimal divergence: a green Argon laser beam starting with a diameter of one centimeter extends to a diameter of three centimeters over a 500-meter journey.





The laser light propagates indefinitely in a well-defined direction, unlike the light of a regular incandescent bulb that emits light in all directions.



Properties of laser light: Monochromatic

Laser radiation always has the same frequency/wavelength while an incandescent bulb emits radiation consisting of photons with different energy.





Monochromatic is linked to the generation of photons by decaying electrons between defined levels and hence with fixed and constant energy. Different active mediums/dopants generate different wavelengths (193 nm - 12000nm)



Laser light properties: Spatial and temporal coherence

Spatial coherence is a consequence of the structure of the resonant cavity of the laser made of two parallel mirrors located at a distance such as to keep the traveling beam in phase with the photons extracted during the stimulated emission process. Time coherence results from the simultaneous extraction of all the photons during the stimulated emission process.



Emissione incoerente: i fotoni vengono emessi casualmente, in tempi diversi e con fase diverse mm annan MAAAA wa mm Emissione coerente: i fotoni vengono emessi simultaneamente e con la stessa fase Coerenza spaziale: le onde hanno la stessa fase in tutti punti della sezione del fascio

> *Coerenza temporale:* le onde conservano la stessa fase nel tempo





Laser Light Properties: Brightness

An important applicative property of laser and consequent of the first two is the defined as the power emitted by the surface unit under a solid unit angle Ω

laser are <u>coherent, monochromatic and with very high brightness</u> light sources

Overview of Laser models: medium

	Laser Name	Oscillation Wavelength	Application
Gas Laser	Helium-Neon Laser	Red/Monochrome 632.8 nm	Optical Axis Alignment Adjusted Length Measurement
	Argon Ion Laser	Blue to Green Multi-line	Optical Axis Alignment, Laser Printing, High-speed Camera Light Source
	Carbon Gas Laser	Infrared 10.6 µm	Metal Welding, Fusing, Processing
	Excimer Laser	Ultraviolet 126 nm to 351 nm	Polymer Micro-processing, Light Source for Academic Purposes (LIF)
	Nitrogen Laser	Ultraviolet 337 nm	Low-priced UV Laser
Solid-state Laser	Ruby Laser	Red 694.3 nm	Holography
	YAG Laser	Red 1064 nm	Metal Micro-processing, Light Source for Academic Purposes (LIF), High-speed Camera Light Source
	Glass Laser	Infrared 1.06 to 1.08 µm	Holography
	Nd (Neodymium) Laser (Nd: YAG, YLF, YVO4, YAIO3)	1064 nm 1047 nm 1053 nm	Optical Axis Alignment / Laser Excitation Micro-processing / Stage Display Light Source
	Titanium: Sapphire	660 nm to 1180 nm	Variable Wavelength Laser
	Fiber Laser	1050 nm to 1620 nm	Long Distance Communication High Temperature Processing
Metal Laser	Helium-Cadmium Laser	Blue White	Medicine Laser Printing
	Copper Vapor Laser	Dual-wavelength 511, 578 nm	Stroboscopic Light Source for High- speed Cameras Uranium Enrichment Pump Laser, Metal Micro-processing
	Gold Vapor Laser	Red	For Medicine, Skin Therapy
Semiconductor Laser	Semiconductor Laser	Red to Infrared	Communication, Solid-state Laser Excitation Light Source, High-speed Camera Light Source, Metal Processing, Laser Pointer Optical Pickup Light Source
Liquid Laser	Dye Laser	300 nm to 1200 nm	Variable Wavelength Laser



Overview of Laser models: pulse duration

Pulsed Lasers

A pulse forming mechanism is needed otherwise lasers run "continuous wave" (CW)

Three types of pulsed operation

- Gain switched (micro or millisecond pulses typically) turn gain on and off (flash lamps, modulate pump)
- Q-switched (nanosecond pulses) modulate cavity loss on times scales > round trip time
- Modelocked (picosecond to femtosecond pulses) modulate cavity loss periodically at roundtrip time



Initial short pulse A pair of gratings disperses the spectrum and stretches the pulse by a factor of a thousand Short-pulse oscillator The pulse is now long and low power, safe for amplification High energy pulse after amplification Power amplifiers Resulting high-energy, ultrashort pulse A second pair of gratings reverses the dispersion of the first pair, and recompresses the pulse.

Simplified scheme for mode-locked fs laser



Fiber laser: a breakthrough in high power laser tech

Unlike most other types of lasers, the laser cavity in fiber lasers is built monolithically by fusion splicing different types of fiber; Fiber Bragg gratings replace conventional mirrors to provide optical feedback. Fiber lasers are pumped by semiconductor laser diodes and can be CW or pulsed. They are by far more compact, electrically efficient in respect to Nd:YAG technology









Femtosecond Lasers for Industry and Science

FEATURES

- 190 fs 10 ps tunable pulse duration
- 2 mJ maximum pulse energy
- 20 W output power
- 1 kHz 1 MHz tunable base repetition rate
- Pulse picker for pulse-on-demand operation
- Rugged, industrial grade mechanical design
- Automated harmonics generators (515 nm, 343 nm, 257 nm, 206 nm)



THE TUNABLE LASER LIGHT SOURCE C-WAVE



SPECIFICATIONS IR b) VIS^{a)} 450 -900 -1300 nm^{b)} Wavelength range 650 nm a) ± 1 nm ± 2 nm Wavelength selection computer controlled Accuracy of wavelength setting internal ± 1 nm ± 2 nm · with external wave-length < 1 MHz °) measurement Power · with 1.5 W > 80 mW > 200 mW pump laser · with 5 W > 200 mW > 400 mW pump laser Amplitude noise < 5 %) <1%9 Beam polarization > 1000:1 TEM₀₀, M² < 1.2^d) Beam profile 0.2 mm c) Beam radius (1/e²) 0.5 mm^{c)} Divergence 0.5 mrad °) 2 mrad^{c)} Linewidth <1 MHz *) Mode-hop-free tuning > 20 GHz °) > 10 GHz *)

Laser applications

Research and industry have made extensive use of laser properties, employing them in the most diverse fields and applications. Its peculiar characteristics make it a flexible, selective and accurate tool and, at the same time, extremely powerful.

Unidirectional and collimated:

- Rangefinders
- Alignment
- Pointing
- Atom cooling

Brightness (High Energy Density):

- Cutting, marking and welding
- Fiber coupling (and associated benefits)
- Nuclear fusion

Monochromatic:

- Selective machining
- Excitation of defined quantum states

DPTOPRI

- Spectroscopy
- Atomic clocks

Coherence:

- Telecommunications
- Metrology

(Gravitational Waves Detection)

- Punctual monitoring variation of physical quantities over wide extended areas (fiber sensing)

The new application frontier is

electrons with photons to increase computing speed



Thank you

Monza: tel 039 834977 Roma: tel 06 87657838



www.optoprim.it info@optoprim.it

Salvatore Salerno: salerno@optoprim.it mobile: 349 2963301