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To initiate a working cycle using this interface, the operator must follow the steps below:

- ENABLE** → Enables the current window.
- 1/3 Pellet Mode** → Selects the operating mode: either 1-pellet or 3-pellet mode.
- Start** → Initiates the Homing phase.
- Insert the pellets into the initial loading chute of the IRIS system.
- Pellet Loaded** → Enter the number of pellets inserted into the initial IRIS loading chute.
- Charge Pellet** → Begins the pellet loading phase into the vacuum chamber. Upon completion, the chamber is coupled.
- Vacuum Request and Collimator Positioning** → This request cannot be performed via the GUI.
- Set the irradiation time(s) for the pellet(s) in the "**OPERATION & SETTING**" section.
- Irradiation** → This phase is executed only if all required operational conditions are met. Pellets are moved one at a time to the irradiation position, where they remain for the duration specified by the user. At the end of the irradiation, an automatic request is sent to the MPS to allow decoupling of the chamber. Once decoupled, the discharge procedure starts, and the first pellet is transferred to the offline measurement position.
- Start/Stop Measurement** → The operator must manually start and stop the measurement for each pellet. Data acquisition from the detectors is **not** managed by the GUI. Instead, the user must use **CoMPASS**, which can be launched via the "**Detector DAQ**" option in the left-hand menu. Once the measurement is started using **Start Measurement**, it can be terminated by pressing **Stop Measurement**. This action discharges the current pellet and positions the next one for measurement.
- Data logging** → In the file Excel save the data related to the timing information displayed in the GUI, as well as the data recorded by CoMPASS.
- Reset** → A reset is required to start a new machine cycle.