

Terzo Incontro Nazionale di Fisica Nucleare INFN 2016

14-16 November 2016, Laboratori Nazionali di Frascati

New Scintillator detectors for nuclear physics experiments

A. Giaz

INFN sezione di Milano



Outline


- ✓ Introduction
 - Nuclear physics experiments
 - Scintillator detectors in nuclear physics experiments
- ✓ Characterization measurements on scintillators
 - $\text{SrI}_2:\text{Eu}$
 - CrBr_3
 - $\text{GYGAG}:\text{Ce}$
- ✓ Measurements of neutron spectra with CLYC scintillator
 - Monochromatic neutrons
 - Neutrons from $^{241}\text{Am}/^9\text{Be}$ source
- ✓ Summary: advantages and disadvantages of each scintillator

Nuclear physics experiments

Nuclear physics experiments, especially with radioactive beams allows to explore unknown mass regions.

A tool to study nuclear structure properties is the gamma decay of GDR (Giant Dipole Resonance).

GDR can be built on excited nucleus (usually fusion-evaporation reaction and compound nucleus) or on ground state.



From the GDR built on excited nuclei:

- ✓ Nuclear shape and deformation
- ✓ Fusion dynamics → dynamical dipole
- ✓ Isospin mixing measurement (CKM-matrix – V_{ud} term)
- ✓ Thermal fluctuations

Scintillators in nuclear physics experiments

Detector requirements:

- ✓ Measurement of high energy gamma rays (~ 15 MeV) \rightarrow Good efficiency
- ✓ Good Time resolution
- ✓ Imaging properties to reduce Doppler Broadening
- ✓ Energy resolution is not mandatory but very useful for:
 - calibration
 - measurement and studies of discrete structures
- ✓ Possibility to discriminate between gamma rays and neutrons using TOF and PSD

Scintillators are the best candidates for this kind of experiments

Material	Light Yield [ph/MeV]	Emission λ_{\max} [nm]	En. Res. at 662 keV [%]	Density [g/cm ²]	Principal decay time [ns]
NaI:Tl	38000	415	6-7	3.7	230
CsI:Tl	52000	540	6-7	4.5	1000
LaBr ₃ :Ce	63000	360	3	5.1	17
SrI ₂ :Eu	80000	480	3-4	4.6	1500
CeBr ₃	45000	370	\sim 4	5.2	17
GYGAG	40000	540	<5	5.8	250
CLYC:Ce	20000	390	4	3.3	1 CVL 50, \sim 1000





The $\text{SrI}_2:\text{Eu}$ scintillator (2'' x 2'')

A. Giaz et al., NIM A 804, (2015), 212

- Energy resolution of $\sim 3.2\%$ at 662 keV
- Slow detector (fall time $\sim 7 \mu\text{s}$)
- Large volume crystals (2'' x 2'') available
- Self absorption

Characterization measurements:



- ✓ Energy resolution up to 9 MeV
- ✓ Crystal scan along the three axes
- ✓ Study of the signal shape





The CeBr₃ scintillator (2'' x 3'')

A. Giaz et al., NIM A 804, (2015), 212

- Energy resolution of $\sim 3.5\%$ at 662 keV
- Very similar to Labr₃:Ce
- Large volume crystals (3'' x 3'') available
- No internal activity



- Characterization measurements:
- ✓ Energy resolution up to 9 MeV
 - ✓ Crystal scan along the three axes
 - ✓ Study of the signal shape




The GYGAG:Ce Scintillator (2'' x 0.7'')

A. Giaz et al., NIM A 804, (2015), 212

- Energy resolution of $\sim 4.0\%$ at 662 keV
- emission wave length 530 nm (yellow)
- No commercially available
- Ceramic detector: any size and shape



Characterization measurements:

- ✓ Energy resolution up to 9 MeV
- ✓ Crystal scan along the three axes
- ✓ Study of the signal shape



Photocathode optimized for the blue region.




A sensor optimized for yellow light will provide a better resolution (SiPM).



The CLYC scintillator ($\text{Cs}_2\text{LiYCl}_6:\text{Ce}^{3+}$)

The elpasolite crystals (CLYC, CLLC, CLLB, TLYC and TLLB) were developed approximately 10 years ago. CLYC:Ce ($\text{Cs}_2\text{LiYCl}_6:\text{Ce}$) belongs to this class of scintillators.

- ✓ Density of 3.3 g/cm^3 , light yield of 20 ph/keV, high linearity, especially at low energy.
- ✓ Energy resolution at 622 keV $< 5\%$ and time resolution of 1.5 ns.
- ✓ Excellent performances in terms of gamma and neutron detection.



Neutron detection

Fast neutrons:

- ✓ $^{35}\text{Cl}(n,p)^{35}\text{S} \rightarrow$ Q-value = 0.6 MeV $\sigma \approx 0.2$ barns at $E_n = 3$ MeV
- ✓ $^{35}\text{Cl}(n,\alpha)^{32}\text{P} \rightarrow$ Q-value = 0.9 MeV $\sigma \approx 0.01$ barns at $E_n = 3$ MeV

$E_{p/\alpha} = (E_n + Q) q_{p/\alpha} \rightarrow$ p or α energy is linearly related to n energy \rightarrow CLYC is a neutron spectrometer

$E_n > 6$ MeV other reaction channels on detectors isotopes \rightarrow not easy neutron spectroscopy

Thermal neutrons:


- ✓ $^6\text{Li}(n,\alpha)t \rightarrow$ Q-value = 4.78 MeV $\sigma = 940$ barns at $E_n = 0.025$ eV.

To fast neutron detection:

^7Li ($^7\text{Li} > 99\%$) enriched CLYC \rightarrow CLYC-7

The kinetic energy of the neutrons can be measured via:

- 1) Time of Flight (TOF) techniques.
- 2) The energy signal



To fast neutron detection:

^6Li ($^6\text{Li} = 95\%$) enriched CLYC \rightarrow CLYC-6




Two measurements:

- ✓ Monochromatic neutrons
- ✓ Continuous neutron spectrum of an $^{241}\text{Am}/^9\text{Be}$ source

Fast Neutron Detection with CLYC

A. Giaz et al., NIM A 825, (2016), 51




Proton Energy [MeV]	Detector Angle	Neutron Energy [MeV]
5.5	0°	3.83
5	0°	3.33
4.5	0°	2.83
5.5	90°	2.68
5	90°	2.30
4.5	90°	1.93



Fast Neutron Detection with CLYC

A. Giaz et al., NIM A 825, (2016), 51

Proton Energy [MeV]	Detector Angle	Neutron Energy [MeV]
5.5	0°	3.83
5	0°	3.33
4.5	0°	2.83
5.5	90°	2.68
5	90°	2.30
4.5	90°	1.93





Continuous neutron spectra

To simulate a continuous neutron spectra the time vs energy matrices (gated on PSD) of the different measured energies were summed.

The pictures show that is possible to separate the contribution of the two different reactions.


The red region includes both the contribution of $^{35}\text{Cl}(n,\alpha)^{32}\text{P}$ reaction and both $^{35}\text{Cl}(n,p)^{35}\text{S}^*$ reaction (1st excited state energy = 1572 keV)



A. Giaz et al., NIM 825, (2016), 51

$^{241}\text{Am}/^9\text{Be}$ Source

$^{241}\text{Am}/^9\text{Be}$ source:




^{12}C can be in different states:

✓ Ground state : $Q = 5.7 \text{ MeV}$

✓ **1st excited state: $Q = 1.3 \text{ MeV}$, $E_\gamma = 4.439 \text{ MeV}$**

✓ 2nd excited state: $E_{\text{th}} = 2.8 \text{ MeV}$ $E_\gamma = 7.654 \text{ MeV}$


✓ 3rd excited state: $E_{\text{th}} = 5.7 \text{ MeV}$ $E_\gamma = 9.641 \text{ MeV}$







Neutron spectra measured in coincidence with a 4.439 MeV γ ray using the TOF technique.



J. Scherzinger, et al., Appl. Rad. and Isotopes, 98, (2015), 74



Measurement of the $^{241}\text{Am}/^9\text{Be}$ spectrum



PDS to separate neutrons from gammas.

$$E_n = E_{\text{mis}}/q - Q$$

$E_n < 7$ MeV: dominant reaction is $^{35}\text{Cl}(n,p)^{35}\text{S}$
 till $E_n < 4$ MeV, for higher energies it is
 necessary to separate different
 contributions. \rightarrow using TOF techniques.

Summary

- $\text{LaBr}_3:\text{Ce}$
 - ✓ Most used, large volume detector available (3.5'' x 8''), fast scintillator, good energy and time resolution, good efficiency
 - X Huge internal activity
- $\text{SrI}_2:\text{Eu}$
 - ✓ Large volume detector available (2'' x 2''), good energy resolution
 - X Self-absorption → better to have small crystals, slow scintillator → no high rate experiment
- CeBr
 - ✓ Similar to $\text{LaBr}_3:\text{Ce}$, fast scintillator, good energy resolution and time resolution, no internal activity
 - X Energy resolution slightly worse than $\text{LaBr}_3:\text{Ce}$
- $\text{GYGAG}:\text{Ce}$
 - ✓ Ceramic scintillator, possibility to built it in every shape and dimension, good energy resolution
 - X Samples not yet available
- $\text{CLYC}:\text{Ce}$
 - ✓ Large volume detector available (3'' x 3''), good energy and time resolution, possibility to discriminate and measure gamma rays and neutrons
 - X Low density and low light yield, two component of signals (long tails)

Thank you for the attention!

Slide di back up

The LaBr₃:Ce scintillator


- Energy resolution of 2.7% at 662 keV
- Fast detector (time resolution ~ 300 ps)
- Good efficiency $\rho = 5.1 \text{ g/cm}^3$
- Good linearity at low energy
- Large volume crystals (3.5'' x 8'') available
- Internal activity



Characterization work already done →
Arrays of LaBr₃:Ce

Example of arrays with LaBr₃:Ce :

- HECTOR⁺
- OSCAR
- PARIS




Position sensitivity with LaBr₃:Ce




Shield PMT with
collimated source

Position of
the centroid




Position reconstruction on event by event



Distance 4 cm



Distance 3 cm



Distance 2 cm






A. Giaz et al., NIM A 772, (2015), 103 & N. Blasi et al., NIM A 839 (2016), 23.

Fast Neutron Detection with CLYC

A. Giaz et al., NIM A 825, (2016), 51




Proton Energy [MeV]	Detector Angle	Neutron Energy [MeV]
5.5	0°	3.83
5	0°	3.33
4.5	0°	2.83
5.5	90°	2.68
5	90°	2.30
4.5	90°	1.93



Fast Neutron Detection with CLYC

A. Giaz et al., NIM A 825, (2016), 51




Proton Energy [MeV]	Detector Angle	Neutron Energy [MeV]
5.5	0°	3.83
5	0°	3.33
4.5	0°	2.83
5.5	90°	2.68
5	90°	2.30
4.5	90°	1.93



Fast Neutron Detection with CLYC

A. Giaz et al., NIM A 825, (2016), 51




Proton Energy [MeV]	Detector Angle	Neutron Energy [MeV]
5.5	0°	3.83
5	0°	3.33
4.5	0°	2.83
5.5	90°	2.68
5	90°	2.30
4.5	90°	1.93



Fast Neutron Detection with CLYC

A. Giaz et al., NIM A 825, (2016), 51

Proton Energy [MeV]	Detector Angle	Neutron Energy [MeV]
5.5	0°	3.83
5	0°	3.33
4.5	0°	2.83
5.5	90°	2.68
5	90°	2.30
4.5	90°	1.93



CeBr₃ – LaBr₃:Ce – NaI:Tl (3'' x 3'')

