


Recent developments on electromagnetic moment measurements using plunger techniques at ALTO and LNL


Alain Goasduff

INFN, Laboratori Nazionali di Legnaro, Legnaro, Italy

Joint LIA COLL-AGAIN COPIGAL POLITA WORKSHOP
April 27, 2016



25 HPGe (GASP) + AC + Complementary detectors



See D. Megoni's talk on Thursday afternoon

GALILEO complementary detectors

Study of weak reaction channels with stable beams :
High efficiency + High resolving power + High CR capabilities

- Light charged particle detectors :
Euclides, TRACE
- Neutron detector :
NeutronWall
- High-energy γ -ray detector :
LaBr₃ detectors
- Heavy-ion detectors :
RFD, Spider¹



Commissioned dets
Near future dets

¹ See K. Hadynska's talk

GALILEO complementary detectors

Study of weak reaction channels with stable beams :
High efficiency + High resolving power + High CR capabilities


- Light charged particle detectors :
Euclides, TRACE
- Neutron detector :
NeutronWall
- High-energy γ -ray detector :
LaBr₃ detectors
- Heavy-ion detectors :
RFD, Spider¹
- Lifetime measurements :
New dedicated plunger IKP-LNL




Commissioned dets
Near future dets

¹ See K. Hadynska's talk


Why measuring electromagnetic moments (EM)?




Why measuring electromagnetic moments (EM)?




Why measuring electromagnetic moments (EM)?




Why measuring electromagnetic moments (EM)?



How to measure them?




How to measure them?



The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,
- Varying the distances between target and stopper,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,
- Varying the distances between target and stopper,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,
- Varying the distances between target and stopper,
- Positioning and alignment of the plunger device,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,
- Varying the distances between target and stopper,
- Positioning and alignment of the plunger device,
- Coupling to complementary detectors,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,
- Varying the distances between target and stopper,
- Positioning and alignment of the plunger device,
- Coupling to complementary detectors,
- The γ -transparency,




Collaboration IKP Cologne - INFN LNL - Based on OUPS developments

The plunger device : Building a new plunger for GALILEO

Essential ingredients :

- One Target and one Stopper/Degrader,
- Stretching both target and stopper foils,
- Varying the distances between target and stopper,
- Positioning and alignment of the plunger device,
- Coupling to complementary detectors,
- The γ -transparency,
- Feedback-system





Collaboration IKP Cologne - INFN LNL - Based on OUPS developments




Commissioning of LNL plunger

Plunger tested with GALILEO and NeutronWall : $^{32}\text{S}(^{154}\text{Sm},6n)^{180}\text{Pt}$




- ^{154}Sm target (1 mg/cm²) on 2 mg/cm² Ta
- Au stopper (10 mg/cm²)



First gated spectra of ^{180}Pt




First gated spectra of ^{180}Pt



Static EM measurements : RIB geometry of the TDRIV technique

A. E. Stuchbery et al. Phys. Rev. C **71**, 047302 (2005)





First test experiment at ALTO :

- **Beam** : ^{24}Mg @ 120 MeV, 0.3 pA
- **Target** : 2.4 mg/cm^2 ^{93}Nb
- **Reset foil** : 1.7 mg/cm^2 Au
- ORGAM γ -ray array : 13 detectors
- 8-fold segmented plastic detector
- Plunger OUPS


Static EM measurements : RIB geometry of the TDRIV technique

A. E. Stuchbery et al. Phys. Rev. C **71**, 047302 (2005)




Angular correlation analysis

$$W(\theta_p, \theta_\gamma, \Delta\Phi, t) = \sum_{kq} a_{kq}(\theta_p) \mathbf{G}_k(\mathbf{t}) D_{q0}^{k*}(\Delta\Phi, \theta_\gamma, 0)$$



Courtesy of A. Stuchbery and A. Kusoglu


Angular correlation analysis



$$R(t) = \left(\prod_{i=1}^n \frac{W_{i\uparrow}}{W_{i\downarrow}} \right)^{1/n}$$



Courtesy of A. Stuchbery and A. Kusoglu

Results from ALTO and the OUPS Plunger



A. Kusoglu et al. PRL **114** 062501 (2015)

Results from ALTO and the OUPS Plunger



A. Kusoglu et al. PRL **114** 062501 (2015)

Extension to Na-like electronic configuration


Method works for H-like ion




Does it work for ions with more than 1 electron ?

Proposed experiment with Na-like ions using $^{56}\text{Fe}^{15+}$:



- Is it possible to observe oscillations ?
- (if yes) Precise and independent g-factor measurement of $^{56}\text{Fe}^{2+}$
- (if yes) Confront ab initio hyperfine field calculations







- Development of a ^{12}C for plunger
 - $230\mu\text{g}/\text{cm}^2$ ^{12}C on $0.5\mu\text{m}$ Ni
- Beam energy tuned to optimize the charge state distribution
- Several conditions to clean the spectrum from :
 - Reactions with the C foil
 - Coulomb excitation in the Ni foil

Preliminary results



- Development of a ^{12}C for plunger
 - $230\mu\text{g}/\text{cm}^2$ ^{12}C on $0.5\mu\text{m}$ Ni
- Beam energy tuned to optimize the charge state distribution
- Several conditions to clean the spectrum from :
 - Reactions with the C foil
 - Coulomb excitation in the Ni foil
- Angular correlations for ORGAM detectors :
 - For stopped and in-flight components
 - For all plunger distances

Preliminary results



- Development of a ^{12}C for plunger
 - $230\mu\text{g}/\text{cm}^2$ ^{12}C on $0.5\mu\text{m}$ Ni
- Beam energy tuned to optimize the charge state distribution
- Several conditions to clean the spectrum from :
 - Reactions with the C foil
 - Coulomb excitation in the Ni foil
- Angular correlations for ORGAM detectors :
 - For stopped and in-flight components
 - For all plunger distances
- Analysis of MINIBALL performed segment by segment

- **OUPS plunger** has demonstrate its reliability during the ORGAM / MINORCA campaigns at ALTO and AGATA@GANIL
- We have successfully build and commissioned a dedicated **GALILEO plunger device** :
 - possibility to accommodate up to 30 Si telescopes from EUCLIDES array
 - possible coupling with RFD for differential plunger technique
- First experimental campaign with this new plunger is on-going
- Development of new forward plastic detector will open the way to high precision g-factor measurement at LNL taking advantage of the ALPI accelerator

THANK YOU VERY MUCH