



Real-Time Quality Assurance in Particle Therapy: simulation and DAQ challenges

Francesco Pennazio
Università di Torino and INFN, TORINO
francesco.pennazio@unito.it

Outline

- Monitoring in HT
- The INSIDE experiment
- INSIDE DAQ system
- INSIDE data processing
- DAQ upgrades
- Monte Carlo Simulations for in-beam PET
 - treatment simulation
 - Detector simulation
- Simulation speed-up

Hadrontherapy monitoring


Why HT monitoring?

Uncertainty in particle range:

- validate the conformity actual dose delivery wrt treatment planning
 - Ensure treatment effectiveness
 - Healthy tissues sparing


HT monitoring (possibly real-time) effects:

- More flexibility in HT treatment planning nearby critical organs
- Higher dose/fraction (hypofractionated plans):
 - More patients per day treated
 - Significant improvement in patient comfort




Hadrontherapy monitoring

Real time monitor of Bragg Peak depth with passive signals by secondary particles




β^+ activity distribution

prompt secondary particles emission




J Pawelke et al., Proceedings IBIBAM, 26.-29.09.2007, Heidelberg




L. Piersanti et al. Phys. Med. Biol. 59 1857


Secondary particles in HT and detectors



Monitoring in HT: Operating principle



Monitoring in HT: in-beam PET



PMMA phantom

2 slices $2.7 \times 2.7 \text{ cm}^2$ ($E_1 = 77 \text{ MeV}$, $E_2 = 105 \text{ MeV}$)

Acquisition time = 254 s, absolute comparison (no normalization applied)



Monitoring in HT: prompt photons

160 MeV proton beam

Energy spectrum of prompt “photons”

(J.Smeets et al., IBA)

Courtesy
Alfredo Ferrari



The *InSide* Project


β^+ activity
distribution
**IN-BEAM PET
HEADS**

Prompt secondary
particles emission
**Tracker +
Calorimeter =
DOSE PROFILER**

INnovative **S**olutions for
In-beam **D**osim**E**try in
Hadrontherapy

Designed to:


- ❑ be integrated in the gantry
- ❑ be operated in-beam
- ❑ **provide an IMMEDIATE feedback on the particle range**



InSide @ CNAO

The *InSide* Project

β^+ activity
distribution
**IN-BEAM PET
HEADS**





INnovative **S**olutions for In-beam **D**osim**E**try in Hadrontherapy


Designed to:

- ❑ be integrated in the gantry
- ❑ be operated in-beam
- ❑ **provide an IMMEDIATE feedback on the particle range**

In-beam PET heads





10x 25 x 5 cm³
Distance from the
isocenter=25 cm



256 LFS pixel crystals (3x3x20mm³) coupled one to one
to MPPCs (Multi Pixel Photon Counters, SiPMs).

20xFE board
(4 Tofpet ASIC^[1]
each)



20xFPGA Xilinx
Spartan6-
SP605

[1] Work partly funded by the European Union 7th Framework Program (FP7/ 2007-2013) under Grant Agreement No. 256984 EndoTOFPET-US and supported by a Marie Curie Early Initial Training Network Fellowship of the European Union 7th Framework Program (PITN-GA-2011-289355-PicoSEC-MCNet).

Inside Preliminary results






Image and profile reconstruction

Preliminary – 25/02/2016
NO AFTER TREATMENT


A



B




D




— Data

— Simulation


beam



PMMA phantom
2 slices $2.7 \times 2.7 \text{ cm}^2$ ($E_1 = 77 \text{ MeV}$,
 $E_2 = 105 \text{ MeV}$)
Acquisition time = 254 s




PMMA phantom with hole
2 slices $2.7 \times 2.7 \text{ cm}^2$ ($E_1 = 77 \text{ MeV}$,
 $E_2 = 105 \text{ MeV}$)
Acquisition time = 155 s



PMMA phantom
Real treatment plan ($E_{\min} = 76.6 \text{ MeV}$,
 $E_{\max} = 124.5 \text{ MeV}$)
Acquisition time = 133 s

The INSIDE in-beam PET DAQ chain



Data processing

Tasks to performed online:

- Data quality monitoring (Gui)
- Energy threshold and **coincidences finding** (high performance dedicated machine and code)
- In-spill and inter-spill discrimination (if applicable)
- MLEM 3D image reconstruction – as soon as enough data is collected
- On-line comparison with expected image

DAQ development and upgrades

Software:

- Developed since project start
- Common C++ framework for analysis of simulation and data
 - > flexible and OO
- Crucial for system check and debug
- Feedback and new functionalities found and implemented on firmware


Firmware:

- Developed since project start
- Basic functions implemented first
- New functionalities added later (i.e. data decoding, calibration and energy threshold)
 - > lower data rate to be processed online by the coincidence finder


Future possible improvements:

- Next-gen (7) FPGA
- System-on-chip Rx layer
 - Fine TDC calibration LUT applicable (consequences: TOF)
- On-line accelerator information brought to DAQ -> finer data matching

In-beam PET simulations



In-beam PET heads simulated geometry



- Simulation time: about 8 h (32-core)
- It is all about number of CPUs (low memory/disk requirements)

In-beam PET simulations: future challenges

Goal: improve **simulation speed**

- Treatment is based on pencil-beams: **easy parallelization**
- Simulation simplification: **bias**, more aggressive energy cuts, simplified transport
- Use of new non-standard tools (planIT?)
- GPU for image reconstruction

Goal: improve **simulation accuracy**

- More accurate reproduction of beam structure (i.e. beam delivery time may not be exactly known a priori)
- Accurate bias and simulation of significant processes other than beta+ decay



Ultimate Goal: fast detector response simulation

- On-line accurate monitoring for cyclotrons and fast duty cycle synchrotrons
- Easier deconvolution of data from background during beam delivery

Conclusions

- Real-time monitoring is a key to fully exploit hadrontherapy
- Different techniques with strengths and pitfalls
- All need reference data to be compared with
- Good preliminary indications from INSIDE in-beam PET test with phantoms
 - Data acquisition in-beam without modification of CNAO clinical workflow
 - Accurate range measurement
- Data acquisition and process require (ongoing) effort on features design
- Simulation validation and accuracy is crucial for HT monitoring
- Simulation speed is a key to extent its purposes

Acknowledgements

Giuseppe Battistoni
Vincenzo Patera
Valeria Rosso

INFN-Torino

Simona Giordanengo
Stefano Bagnasco

INSIDE PET team

Maria Giuseppina Bisogni
Niccolo' Camarlinghi
Alberto Del Guerra
Veronica Ferrero
Elisa Fiorina
Giuseppe Giraudo
Matteo Morrocchi
Francesco Pennazio
Cristiana Peroni
Maria Antonietta Piliero
Giovanni Pirrone
Richard Wheadon

ASIC designers

Angelo Rivetti
Manuel Dionisio Rolo

CNAO personnel

Marco Donetti
Mario Ciocca
Marco Pullia

INFN – Torino Technical Staff

Giuseppe Alfarone
Florea Dumitrache
Marco Mignone
Barbara Pini
Mauro Scalise
Antonio Zampieri