From MACACO to MACACO II: a Compton telescope for hadron therapy treatment monitoring.

Gabriela Llosá, John Barrio, Ane Etxebeste, Laura Granado, Carlos Lacasta, Enrique Muñoz, Josep F. Oliver, Ana Ros, Carles Solaz.

Instituto de Física Corpuscular IFIC (CSIC-UV)/ IFIMED

IRIS Group: Image Reconstruction, Instrumentation and Simulations for Medical Imaging

http://ific.uv.es/iris

Gabriela.llosa@ific.uv.es



Gabriela Llosá

Outline

- Introduction to hadron / particle / ion beam therapy.
- Hadron therapy monitoring.
- Compton camera.
- Compton telescope @ IFIC/IFIMED
 - Final results with MACACO
 - First results with MACACO II

MACACO: Medical Applications CompAct COmpton camera



Hadron therapy

• Charged particles (protons or carbon ions) instead of photons.



Gabriela Llosá

Advantages

- Precise delivery to tumour area => increase of cure rates and reduction of side and long term effects and secondary cancer.
- Higher relative biological effectiveness (RBE) than photons.
- Large benefit over conventional radiation therapies in some cases (ocular tumours, children, organs at risk, radioresistant tumours).



Gabriela Llosá

Hadron Therapy in the world

- Gaining increasing importance -> growing number of centers in the world.
- >70 centers in operation; > 40 under construction; ~20 planned.



Hadron therapy centers in the world (from http://ptcog.web.psi.ch/ptcentres.html)

Gabriela Llosá

But...

• High precision is not an advantage if the aim is not correct (machine gun vs. precision rifle).



- Large margins have to be applied in treatment planning.
- Range uncertainties are considered one of the main obstacles for a wider application of this techniques

Real-time monitoring is essential

Gabriela Llosá MEDAMI 2017, Orosei, 30 May - 4 June

ION BEAM

DOSE DEPOSITION BY ELECTRONS THROUGH IONIZATION (electromagnetic interaction)

EMISSION OF SECONDARY PARTICLES (nuclear interaction)

IRRADIATED TISSUE ATOMS







- Different mechanisms, but correlated
- Indirect estimation of dose deposition- can be used for monitoring

Gabriela Llosá



Gabriela Llosá MEDAMI 2017, Orosei, 30 May - 4 June



Gabriela Llosá

Current monitoring - PET

- Positron Emission Tomography (PET) + MC currently employed.
- Verification: comparison of β + activity predicted from planned dose and measured with PET.





Gabriela Llosá

PET limitations

- PET Limitations:
 - Positron production does not follow irradiation immediately.
 - Biological washout- activity carried away by metabolic processes.
 - Low amount of β + activity induced- low efficiency.
 - Difficult online studies partial ring.
 - Photons produce significant background.
- PET improvements:
 - Improved systems.
 - Models for washout.
 - Use of short-lived isotopes.
 - TOF PET to minimize gap effects.



Gabriela Llosá



Gabriela Llosá

Prompt gammas

- Prompt gammas also emitted from nuclei excited during therapy and correlated with dose.
 - MeV in a continuous energy spectrum. High energies.
 - Emitted within ns after irradiation.
 - Higher amount than positron emitters.



Gabriela Llosá

How to detect them?



F. Roellinghoff et al. Phys. Med. Biol. 59 (2014)

J. Smeets. Phys. Med. Biol. 57 (2012)

Gabriela Llosá

Compton Imaging



The cone surface is projected to the reconstruction volume.

The intersection of several cone surfaces yields the position of the source.

Backprojection





$$\lambda_j^{n+1} = \frac{\lambda_j^n}{s_j} \sum_{i=0}^N \frac{t_{ij}}{\sum_k t_{ik} \lambda_k^n}$$

Gabriela Llosá MEDAMI 2017, Orosei, 30 May - 4 June

Compton camera configuration



Detector 2



Detector 3

Gabriela Llosá

Scatterer + absorber: 2 interactions

Problems if the photon energy is unknown or if it can escape (MeV)

$$cos\theta = 1 - m_0 c^2 \left(\frac{1}{E_0 - E_e} - \frac{1}{E_0}\right)$$

Multilayer: 3 interactions in 3 detectors (+ correct ordering):

- Energy determined
- lower efficiency

$$\cos(\theta) = 1 - \frac{E_1 m_e c^2}{E_0 (E_0 - E_1)}$$

$$E_0 = E_1 + \frac{1}{2} \left(E_2 + \sqrt{E_2^2 + 4 \frac{E_2 m_e c^2}{1 - \cos \theta_2}} \right)$$

Compton telescope @ IRIS



- Combination of :
 - 2 int events (high efficiency) +
 - 3 int events (high resolution).

+ Estimate both position and energy simultaneously from 2interaction events

NO ABSORPTION REQUIRED

Gabriela Llosá

Compton telescope @ IRIS



Gabriela Llosá

Detectors

Detector 1



Detectors 2 and 3

Monolithic crystals: E resol: 7-10% FWHM Sp resol: 1.3 mm FWHM Expected DOI:~ 2 mm FWHM,

Gabriela Llosá

DAQ system

- VATA64HDR16 ASIC from IDEAS 64 channels
- Connected to a DAQ system made at IFIC.
- Compact and portable system



Two-layer simulations vs. data



Two planes: different energies



RESOLUTION IMPROVEMENT WITH ENERGYGabriela LlosáMEDAMI 2017, Orosei, 30 May - 4 June

Two planes: different positions

- Sources placed in different positions, separated 20 mm.
- Points reconstructed separately.

GOOD POSITIONING ACCURACY

Na-22 selecting 1275 keV peak



Y-88 selecting 1836 keV peak

Det 2

Det 1

y source⊾

40 mm

80 mm



Gabriela Llosá

Two planes: different configurations



CONF 1 CONF 2 CONF 3 GEOMETRY DETERMINES RESOLUTION TRADE-OFF WITH EFFICIENCY



Na-22 selecting 1275 keV peak MEDAMI 2017, Orosei, 30 May - 4 June

Gabriela Llosá

24

Two planes: Na-22 + Y-88





- Sources separated 40 mm.
- Data from both sources reconstructed together





Gabriela Llosá

Three planes: Na-22

Na-22: 650 keV< Esum < 1350 keV.







Gabriela Llosá

Three planes: Y-88

Y-88: 1100 keV< Esum < 1950 keV.



Gabriela Llosá

Efficiency



Tests in a proton beam

- Tests at KVI-CART, AGOR cyclotron (Groningen).
- Proton beam, 150 MeV, ~ 10^8 prot/sec.
- Graphite and PMMA targets.



Gabriela Llosá

Tests in a proton beam

- Data with two layers in coincidence.
- PMMA target shifted to simulate Bragg peak variations.

SHIFT OBSERVED IN BRAGG PEAK POSITION



Tests with high energy gammas

- Tests at HZDR 3MV Tandetron (Dresden).
- Gammas 4.439 MeV. ${}^{15}N(p,\alpha\gamma){}^{12}C$ reaction.
- Low gamma emission rate.



Gabriela Llosá

Tests with high energy gammas



From MACACO to MACACO II

- System limitations: Energy resolution, timing resolution, electronics low dynamic range.
- New version:
 - New detectors to improve energy resolution.
 - Acquisition of two and three interaction events simultaneously with a new coincidences board.
 - New electronics?
 - Det 3 under development, 4 times larger

New photodetectors:

- TSV, 8x8 pixels.
- Better uniformity.
- Better PDE

Hamamatsu S13361-3050AE-08



Gabriela Llosá

Second version of the system



LaBr3 crystal: Cs137 spectrum







Linearity and energy resolution







Gabriela Llosá

Tests at CNA

- Spanish National Accelerator Centre, Sevilla.
- 18 MeV protons on a graphite target \rightarrow 4.4 MeV gamma rays.



Gabriela Llosá

Conclusions

- Final results obtained with MACACO I. Critical points identified.
- MACACO II is being assembled. First tests carried out.
 - Improved energy resolution.
 - Very promising first results.
 - Other improvements underway → Need to address application requirements.
- Need to carry out tests in clinical beams.

Acknowledgements

- This work is supported in part by:
 - ENVISION project, co-funded by the European Commission under FP7 Grant Agreement num 241851 and ENTERVISION project.
 - Generalitat Valenciana (GV/2013/133).
 - Ministerio de Economía y Competitividad (FPA2014-53599-R).
 - IFIC- Centro de Excelencia Severo Ochoa (SEV-2014-0398).
- Group members are supported by Ramón y Cajal IFIC-Severo Ochoa, UVEG Atracció de Talent and Generalitat Valenciana contracts.

Thank you! Questions?

