

# Status of the HZB Cyclotron

A. Denker, J. Bundesmann, T. Damerow, T. Fanselow,
W. Hahn, D. Hildebrand, U. Hiller, I. Kailouh, S. Seidel,
U. Müller, C. Rethfeldt, J. Röhrich, C. Zimmer
Helmholtz-Zentrum Berlin, Protons for Therapy

- k=130 cyclotron (former VICKSI, ISL)
- two injectors:
  - 2 MV Tandetron<sup>TM</sup>, standard for therapy
  - 6 MV Van-de-Graaff, backup, time structures





- k 130 cyclotron (former VICKSI, ISL)
- two injectors:
  - 2 MV Tandetron<sup>TM</sup>, standard for therapy
  - 6 MV Van-de-Graaff, backup, time structures
- three target stations:
  - treatment room
  - experimental station
     (I<sub>max</sub> = 10 nA)
  - beam line end for tests in cyclotron vault







# • since 2007:

accelerator financed by external means (hospital, third party funds)

- since 2009: accelerator R & D again part of HZB's research portfolio
  - bachelor and master theses, PhD students
  - small budget for research



#### therapy:

- 68 MeV protons, quasi-DC, broad beam (ø 50 mm), I<sub>patient</sub> < 3 nA</li>
- deliverable by either Van-de-Graaff or Tandetron as injector
- experiments:
  - broad or focused beam
  - quasi DC to single pulses with t < 1 ns (single turn extraction)</li>
  - changes in intensity: 0.1 pA  $\leq$  I<sub>target</sub>  $\leq$  1500 nA
  - 68 MeV protons, <sup>4</sup>He: 50 MeV, 75 MeV, 90 MeV
  - <sup>3</sup>He: 50 MeV under development

87% 13%

Therapy
Acc.Dev.
Dosimetry
Med.Phys.
PIXE
Radhard
Nucl.Phys.





- less scheduled beam time: major events → huge impact on statistics
  - in 2015: human error increase of injector voltage too fast
  - many errors appear during start-up of accelerator
  - past years: less than 5% when patients present







- most years: cyclotron, especially RF, is main "culprit"
  - modernize RF electronics: replacement of old low level control with modern system from iThemba labs (poster)
- external power failures
- since 2011: Tandetron = standard injector for therapy







#### until 2006: 14 GHz ECR plus RFQ as injector for heavy ions







# 2007/08: replacement of RFQ with 2 MV tandetron





- mechanical constraints in positioning
  - emergency exit
  - access to cyclotron
- at the end: compromise

# Beam Transport RFQ → Cyclotron





# Beam Transport Tandetron → Cyclotron





# Tandetron → Cyclotron: Tuning Issues

HZB Helmholtz Zentrum Berlin

- start parameters from Tandetron not well known
  - parameters from Cadarache
- electrostatic quadrupole:
- triplet with 3 (three!) power supplies: U-, U+, U-
- asymmetric quadrupole



difficult to calculate
 asymmetric beam

# Tandetron → Cyclotron: Tuning Issues



- start parameters from Tandetron not well known
  - parameters from Cadarache
- electrostatic quadrupole:
  - triplet with 3 (three!) power supplies: asymmetric quadrupole
- asymmetric beam
- was observed on beam profile monitor
- interpreted as broad x and narrow y beam,
- slightly off axis in y
  - possible off set in alignment
- "experimental" tuning: good transmission to and through cyclotron
- Iow priority







- $\rightarrow$  tuning ambiguous
- installation of a harp for better reproducibility
- master thesis of M. Burmeister (now at Stryker)
- 48 wires in x and y (broad beam)
- special soldering unit:  $\Delta$  wire = 0.5 mm

USB microscope







- BPM not at focal point + beam emittance defined close to BPM
- $\rightarrow$  tuning ambiguous
- installation of a harp for better reproducibility
- 25 wires in x and y (broad beam)
- mounted on standard movement unit (150 mm travel)





- BPM not at focal point + beam emittance defined close to BPM
- $\rightarrow$  tuning ambiguous
- installation of a harp for better reproducibility
- 25 wires in x and y (broad beam)
- mounted on standard movement unit
- connection via flat cable (Würth Elektronik) and PCB boards





- BPM not at focal point + beam emittance defined close to BPM
- $\rightarrow$  tuning ambiguous
- installation of a harp for better reproducibility
- 25 wires in x and y (broad beam)
- mounted on standard movement unit
- connection via flat cable and PCB boards
- vacuum feed through: PCB board and epoxy
  - after 6 hours: vacuum better than 2 10<sup>-7</sup> mbar
  - leak tested: 1 10<sup>-9</sup> mbar/(l s)
  - mass spectrometer: nothing dangerous for electrostatic quadrupole nearby





- BPM not at focal point + beam emittance defined close to BPM
- $\rightarrow$  tuning ambiguous
- installation of a harp for better reproducibility
- 25 wires in x and y (broad beam)
- mounted on standard movement unit
- connection via flat cable and PCB boards
- vacuum feed through: PCB board and epoxy
- electronics: "harp box" iThemba Labs







- similar profiles for harp and BPM
- beam as calculated, however: beam does not match into cyclotron
- operators can quantify beam width
  better reproducibility







- perfect transmission through cyclotron:
- in X and Y: harp profile identical to BPM!
- why double peak in Y?
- until now: explained as slight misalignment
- two beams ??





HZB Helmholtz Zentrum Berlin

- perfect transmission through cyclotron
- in X and Y: harp profile identical to BPM!
- why double peak in Y?
- nevertheless: beam fully fulfils requirements of therapy
- data for reverse beam line calculations available



channel



12 16 20 24 28 32 36 40 44 48

channel

8

4

### **Tumour Therapy**

- on-going mission: 12 therapy weeks per year,
   ~ 210 patients/year
- in addition to treatment: master students for medical physics
  - Beuth Hochschule für Technik
  - Martin-Luther-Universität Halle

irradiation of mice (one eye) for investigation of radiation retinopathy







## Conclusion

- reliable accelerator operation, past years: uptime usually > 95 %
- new ion species and energies: all light ions
- new installations
  - harp  $\rightarrow$  beam line calculations now verifiable
  - new LLRF control for the cyclotron
- students for research
- 27. Jan. 2017: 3000. patient
   Thank you for your attention!

1998 1999 2000 2001 2002 2003 2004 2005 2006 2007 2008 2009 2010 2011 2012 2013 2014 2015 2016