



"Monitoring chamber simulation and analysis"



I. Margjeka, N. De Filippis, R. Aly, W. Elmetenawee
Università, Politecnico and INFN of Bari


F. Cuna, A. Corvaglia, G. Chiariello, F. Grancagnolo, M. Manta,
A. Miccoli, M. Panareo, G. F. Tassielli
Università del Salento, INFN, Lecce

- Framework structure
- Simulation on the Monitoring chamber
- Analysis on the simulation's results

- Libraries:
 - Garfield
 - ROOT
- Simulation and analysis framework:
 - GAS-SIM: Simulation different gas-mixtures
 - SIM: Simulation of electric field configuration and tracks through the chamber
 - Theoretical calculation of the drift velocity: `calcMeanVdrift`
 - Analysis macros

Simulation on the Monitoring chamber

- The macro for the simulation is: `sim.C`
 - For symmetries we are performing the simulation on one quarter of 2D-scheme of the chamber
 - Dimensions: $-2\text{cm} < x < 0\text{cm}$, $0 < y < 1\text{cm}$



1) Implementation of the chamber's structure:

- Medium (mixture of gas)
- Geometry: the box in itself
- Wires

2) Electric field configuration:


- 1 We vary the V on the guard-wires ($[-800, 0] V$; $25 V$)
- 2 We vary the V on the sense-wires ($[0, 1200] V$; $25 V$)
- 3 The voltage on the field-wires is fixed on $-2000 V$
- 4 We fix $Oz = 0$ and scan the electric field per each Oxy -plane along the Oy -axis from ($[0, 0.9] cm$; $0.1 cm$)

3) Generation of the tracks:


- 1 Avalanche parameters possible to set up like: time and spatial interval, etc ...
- 2 Tracks: start from $y = 1$ and $x = \pm 1$ and in the negative Oy -direction with an angle θ extracted from a uniform distribution $\theta \pm 12^\circ$

Electric field configuration

Electric field configuration
with $V_g = 0V$, $V_s = 0V$




Electric field configuration with
 $V_g = -350V$, $V_s = 925V$




Potential configuration




Potential configuration
with $V_g = 0V$, $V_s = 0V$



Potential configuration
with $V_g = -350V$, $V_s = 925V$



Chamber Physics parameters




.... Part II:
"Analysis Work"

Analysis on the simulation's results

⇒ Macros used for the analysis on the simulation's results are:


- `analisiSum.C`
- `AnalData1.C`
- `calcMeanVdrift.C`
- `plotVdTrend.C`

- The main goal is to determine the optimized value for V_s and V_g to ensure a high and uniform electric field in the two drift cells.
 - From the `sim.C` we obtain the RMS, skewness, electric field variation ΔE and the electric field configuration per each Oxy -plan along the Oy -axis and for different value of V_s and V_g .
- ⇒ An example of the scan of the electric field, ΔE , RMS, skewness average values on the plane corresponding to the first sense wire in analysis.




- We calculate the weights for the RMS, skewness and the variation of the electric field ΔE :
 - We assign a weight value between 0-1 to each bin of the histograms normalizing the Oz -axis.
 - We divide the value of each bin by a reference value.
 - The reference value is chosen as a mean value to relative distributions with a confidence value of 3σ .
 - If the weight value is > 1 , then the given weight value will be just 1.
 - The goal is to append a bigger weight to the minimal variation of the field, RMS and skewness, by doing a complementary histogram, where each bin is weighted as $1 - weight_{previous}$.

⇒ An example of the scan of ΔE , RMS and skewness for the first sense wire.




- We analyze the information of each scan along the O_y -axis, creating 3 histograms related to the RMS, skewness and ΔE , where each bin is filled with the mean values of all the scans.
- Moreover, we create a global histogram, where each bin contains the average value of the weights from the previous histograms.




- top left: the global histogram that contains the average of the other three histograms.
- top right: average of ΔE for all the plane scans.
- bottom left: average of RMS for all the plane scans.
- bottom right: average of Skewness for all the plane scans.

For the two sense wires in analysis, we study the electric field configuration, the number of collected electron and the gain value.




- We assign a weight also to the histograms for the gain and number of the collected electrons in the sense wires.
- The weights are assigned in the same way as in the previous explanation.
- In this case we divide the value of each bin by the number of simulated tracks.
- In case of the gain:
 - The weights per each bin of the gain histogram has been assigned in logarithmic scale, where each bin has divided by a value of 10^5 (reference value for proportional counting region).
 - If the weights are between 1-2, the new weights are obtained as $2 - weight_{previous}$.
 - If the weights are > 2 , we assign the value 0 for those weights.

⇒ The weighted histograms for gain and number of collected electron of the two sense wires in analysis.



analisiSum.C: Part 4

- At last, we build a global histogram, where each bin contains the average value of the weights from the histograms of:
 - Electric field uniformity
 - Gain
 - Collected electrons of each sense wire
- From this global histogram, we can then obtain the optimized voltage values for sense and guard wires:
 - $V_s = 925V$
 - $V_g = -350V$



- The main goal is to determine the drift velocity values as a function of He content or pressure/temperature variations from the double peak Gaussian distribution.

Some remarks:


$$\begin{cases} t_2 = \frac{t_1 + t_3}{2} \mp \frac{2\delta}{v_d} \\ t_3 = \frac{t_2 + t_4}{2} \pm \frac{2\delta}{v_d} \end{cases}$$

$$\rightarrow t_3 - t_2 \rightarrow \Theta = (t_1 + t_3 - 2t_2) - (t_2 + t_4 - 2t_3)$$

$$\Theta = \begin{cases} \text{Left} : \Theta_+ = +\frac{8\delta}{v_d} \\ \text{Right} : \Theta_- = -\frac{8\delta}{v_d} \end{cases}$$


- Θ will have double gaussian distribution.
- The distance between the two peaks is related to drift velocity:

$$v_d = \frac{16\delta}{\Delta\Theta}, \sigma_{v_d} = \sqrt{\left(\frac{16}{\Delta\Theta}\right)^2 \sigma_\delta^2 + \left(\frac{-16\delta}{(\Delta\Theta)^2}\right)^2 \sigma_{\Delta\Theta}^2} \quad (1)$$



- Identification of the drift time t_d on each single wire:
 - Weighted average of the drift time of those electrons with a minimal collected charge equal to 10^5 electrons.
 - Otherwise the information on the wire will be skipped.
 - The weights are evaluated by the ratio of the charge value recorded on the wire over the reference value of 5×10^5 .

Here an example of the double peak gaussian distribution for a mixture 89% He - 11% iC_4H_{10} , at $P = 760\text{torr}$ and at $T = 300\text{K}$.



The drift value extrapolated from this distribution is:

$$v_d = 1.97 \pm 0.03 \text{ cm}/\mu\text{s}$$

- The main goal is to make a comparison between the drift velocity value from the double peak gaussian distribution and its theoretical value.
- It simulates theoretically the drift velocities v_d for different gas mixtures according to two methods:

- $v_{d,1} = f(|\vec{E}|)$


- $v_{d,2} = \frac{D_d}{t_d}$

Both of the v_d -values are taken as a mean value:

$$\bar{v}_{d,1} = \frac{\int v_d \delta E}{\Delta E}$$

$$\bar{v}_{d,2} = \frac{v_{d,2}}{N_{DriftLines}}$$

- It does the plots of the theoretical and experimental drift velocity:



There is still a shift between theoretical value and experimental value of drift velocity.