

## A dedicated MRI/PET System for Radiotherapy Treatment Simulation

*Sunday, 29 May 2022 10:10 (20 minutes)*

The Unity 1.5T MR-linac has clinically been introduced to provide stereotactic precision for moving targets. We are now developing an integrated 1.5T PET/MRI system with identical MRI performance for planning of MR-linac treatments. The system offers 70 cm bore size and includes PET detectors mounted in the gap of a split gradient coil using the Hyperion III platform. Recently, the first ring of PET detectors was installed. We performed extensive interference tests and acquired the first simultaneous images using a brain phantom. Our vision is that the combination of MR-linac with PET/MRI will allow the use of PET information for stereotactic radiation therapy guidance.

**Primary author:** Dr SCHUG, David (RWTH Aachen University, PMI - Physics of Molecular Imaging; Hyperion Hybrid Imaging Systems GmbH)

**Co-authors:** MÜLLER, Florian (Department of Physics of Molecular Imaging Systems, Institute for Experimental Molecular Imaging, RWTH Aachen University); Dr WEISSLER, Bjoern (RWTH Aachen University, PMI - Physics of Molecular Imaging; Hyperion Hybrid Imaging Systems GmbH); RADERMACHER, Harald (RWTH Aachen University); GEGENMANTEL, Eike (Hyperion Hybrid Imaging Systems GmbH); BEIJST, Casper (UMC Utrecht); BRANDERHORST, Woutjan (University Medical Center Utrecht); Prof. LAGENDIJK, J.J.W. (University Medical Center Utrecht); BORGIO, Martino (Futura Composites); VAN DUIN, F. (Philips Best NL); MOLLINK, Jurgen (Philips Best NL); BORST, Govaert; VERHEYEN, Marc (Philips Best NL); KLOMP, Dennis (University Medical Center Utrecht); SCHULZ, Volkmar (RWTH Aachen University, Hyperion Hybrid Imaging Systems, Fraunhofer MEVIS); DE JONG, Hugo (UMC Utrecht)

**Presenter:** Dr SCHUG, David (RWTH Aachen University, PMI - Physics of Molecular Imaging; Hyperion Hybrid Imaging Systems GmbH)

**Session Classification:** PET/MR and SPECT/MR systems and applications

**Track Classification:** PET/MR and SPECT/MR systems and applications