

Dead Time Effects and Image Quality Evaluation at High Activities for the SAFIR Dual Ring Prototype



A closer look at challenges faced by a high rate PET system

PSMR-TBP 2022

28.5 – 01.06.2022

Elba, Italy

Jan Debus, ETH Zurich
for the SAFIR collaboration



Motivation:

- Preclinical PET insert for quantitative kinetic modeling, e.g. ^{15}O -water time-activity curves
 - Fully MR-compatible, simultaneous image acquisition in mere seconds
- **Small Animal (mice and rats) Fast (image times ~ seconds) Insert for MRI (SAFIR)**

Requirements:

- PET insert operating within a Bruker BioSpin 70/30 MRI scanner (7T)
 - Spatial resolution: ~2 mm, Image acquisition times: ~5 s,
 - Target activity of 500 MBq in mice and rats
- **Result:** SAFIR Dual Ring Prototype (DRP), reduced FOV prototype, completed in 2019




Scanner Layout

- 199 mm outer, 114 mm inner diameter
- Carbon fibre structure, fully air cooled
- 12 sectors, each sector RF screened

Crystal Geometry

- 2.1x2.1mm LYSO crystals, 2.2mm pitch
 - Crystals 1:1 coupled to SiPM Matrices
 - SiPM digitization using PETA6-SE ASIC
 - 16 rings, 12 x 15 crystals per ring
- **2880** crystals in total, 35.6mm Axial FOV



Digital readout

Common electronics

Crystals, SiPM arrays,
Readout ASICs



Operating at 500MBq:

- SAFIR DRP characterized according to NEMA NU-4, including NECR @ 500MBq
- NEMA requires evaluation of Image Quality Phantom at 3.7 MBq
- ➔ Of course, Image Quality at higher activities of interest




Measurements:


- NEMA Image Quality Phantom at 100 MBq up to 500 MBq
- Evaluation of Spill-over-Ratios, Recovery Coefficients and Uniformity according to NEMA

Data processing:

- Data calibration & coincidence sorting performed offline (custom software)
- Coincidence timing window: **500 ps**
- Coincidence energy window: **391 keV – 601 keV**
- Recon: **STIR, OSMAPOS, 12 Subsets, 30 Iterations**
- Corrected for randoms, scatters, normalization, attenuation



Recovery Coefficient
(RC)



Spill-Over-Ratio
(SOR)



Uniformity




Measurements:


- NEMA Image Quality Phantom at 100 MBq up to 500 MBq
- Evaluation of Spill-over-Ratios, Recovery Coefficients and Uniformity according to NEMA

Data processing:


- Data calibration & coincidence sorting performed offline (custom software)
- Coincidence timing window: **500 ps**
- Coincidence energy window: **391 keV – 601 keV**
- Recon: **STIR, OSMAPOSL, 12 Subsets, 30 Iterations**
- Corrected for randoms, scatters, normalization, attenuation





Recovery Coefficient
(RC)




Spill-Over-Ratio
(SOR)




Uniformity




Recovery Coefficient (RC)



Spill-Over-Ratio (SOR)




Uniformity



Notes:

- Meas. times chosen to equate 20 minutes @ 3.7 MBq
- Measurement error due to mechanical movement possible



Notes:

- Meas. times chosen to equate 20 minutes @ 3.7 MBq
- Measurement error due to mechanical movement possible

Conclusions:

- Slight variations in all parameters observable
- No significant degradations detected



Operating at 500MBq:

- SAFIR DRP characterized according to NEMA NU-4, including NECR @ 500MBq
- NEMA requires evaluation of Image Quality Phantom at 3.7 MBq
- Of course, Image Quality at higher activities of interest


On the topic of Data Loss:

- Measured count rates at 500MBq lower than simulated
- Dead-Time effects of detector channels suspected to cause data loss
- Investigate Dead-Time effects at higher activities



Detector Dead time:

- Readout channels enter 'busy' period after hit
 - Photons arriving during that time are not registered
 - 'Deadtime' ranges from 450 ns to 2.4 μ s for PETA6
 - Cross-talk, noise & intrinsic Lutetium activity can trigger channels as well
 - At high event rates this can be a significant factor
- Increased loss of singles & coincidences





Test-Trigger Events:

- PETA-6 supports injection of 'dummy' signals for all channels
 - Will only be read out if channel is not 'busy', else lost
- Injection at e.g. 100 Hz, then recording the number received
 - Estimate percentage of real events lost per channel per second
 - Averaging over the entire detector to estimate data loss


Measurement:

- Cylindrical Calibration Phantom (40 mm x 20 mm)
- Activity: F18 in water, 484 Mbq to 17 MBq



Single Event Loss:

- At 484 Mbq, an average of **11.2 %** of single events were lost
- Near linear dependency






Single Event Loss:

- At 484 Mbq, an average of **11.2 %** of single events were lost
- Near linear dependency

Coincidence Loss:


- For each LOR, channel loss rates multiply
- 21% coincidence loss at 500MBq






Coincidence Loss Correction:

- Multiply individual channel loss rates for each LOR
- Use results as correction factor for LOR count rates
- Correction restores linear dependency of count rate and measurement activity






- Evaluation of Image Quality Phantom measurement performed before and after correcting for data loss
- Comparison of parameters at 100 and 500 MBq shows almost no deviations






- Evaluation of Image Quality Phantom measurement performed before and after correcting for data loss
- Comparison of parameters at 100 and 500 MBq shows almost no deviations





- Evaluation of Image Quality Phantom measurement performed before and after correcting for data loss
 - Comparison of parameters at 100 and 500 MBq shows almost no deviations
- **Correction yields important benefits for quantitative PET measurements with varying activities while not degrading Image Quality**





SAFIR-I built & in use

- Dual Ring Prototype was upgraded
- Axial FOV increased by 50%
- First Time-Activity curves measured

SAFIR-II under construction

- Increased Axial FOV of 142.4 mm (~ x4)
- **New ASIC:** PETA6 → PETA8, reduces power consumption by 50%
- Adjustments to crystal reflectors, reduced cross-talk

Continued Dataloss Investigations

- PETA-8, featuring improved readout logic, decreases deadtime
- Measurements with SAFIR-II, impact on In-Vivo Data
- Integration of deadtime into simulations
- Investigation on count-rate performance after deadtime-correction



The SAFIR Collaboration

- Institute for Particle Physics and Astrophysics, ETH Zürich
- Institute for Pharmacology and Toxicology, University of Zürich
- Department of Nuclear Medicine, University of Zürich
- Institute for Biomedical Engineering, ETH Zürich
- Institute of Computer Engineering, University Heidelberg
- University Medical Center, Groningen

This work was supported by the ETH Zurich Foundation through ETH Research Grant ETH-30 14-2, and by an SSSTC Exchange Grant EG 02-03201 and by the Einrichtungskredit of University of Zurich.



SPARE






- Measurement & evaluation of Image Quality Phantom performed before and after correcting for data loss
- Comparison of parameters at 100 and 500 MBq shows almost no difference

Activity	100 MBq		500 Mbq	
	no	yes	no	yes
SOR (Air)	0.157	0.157	0.163	0.164
SOR (Water)	0.179	0.180	0.192	0.193
Uniformity [%]	3.17	3.17	2.89	2.87
RC (1 mm)	0.175	0.176	0.199	0.205
RC (2 mm)	0.315	0.315	0.325	0.334
RC (3 mm)	0.520	0.520	0.503	0.515
RC (4 mm)	0.688	0.688	0.646	0.661
RC (5 mm)	0.867	0.867	0.876	0.896



Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876





Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876

Notes:

- Measurement error due to mechanical movement possible
- Meas. times chosen to equate 20 min @ 3.7 MBq



Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876

Notes:

- Measurement error due to mechanical movement possible
- Meas. times chosen to equate 20 min @ 3.7 MBq



Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876

Notes:

- Measurement error due to mechanical movement possible
- Meas. times chosen to equate 20 min @ 3.7 MBq



Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876

Notes:

- Measurement error due to mechanical movement possible
- Meas. times chosen to equate 20 min @ 3.7 MBq



Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876

Notes:

- Measurement error due to mechanical movement possible
- Meas. times chosen to equate 20 min @ 3.7 MBq



Activity [MBq]	100	200	300	400	500
Meas. Time [s]	48	24	16	12	10
SOR (Air)	0.157	0.163	0.165	0.163	0.163
SOR (Water)	0.179	0.179	0.182	0.184	0.192
Uniformity	3.17%	3.02%	3.03%	3.12%	2.89%
RC (1 mm)	0.175	0.151	0.190	0.187	0.199
RC (2 mm)	0.315	0.289	0.321	0.316	0.325
RC (3 mm)	0.520	0.476	0.506	0.499	0.503
RC (4 mm)	0.688	0.649	0.662	0.645	0.646
RC (5 mm)	0.867	0.856	0.865	0.869	0.876

Conclusions:

- Slight variations in all parameters observable
- No significant degradations detected