



1

The HERD Space Mission



Timing resolution checks

F. Alemanno, M. Fernandez, F. Barbato, I. De Mitri, A. Di Giovanni, M. Di Santo, D. Kyratzis, A. Parenti, L. Silveri, A. Suma, L. Wu Gran Sasso Science Institute (GSSI) & INFN-LNGS

Goal: to assess the intrinsic timing properties of bar+SiPM only



Components

SiPM 1, SiPM 2: Hamamatsu S13360 – 3025CS

Scintillator coupled to SiPM1: Eljen (EJ 200)

Scintillator coupled to SiPM2: Saint-Gobain (BC-404)

Linear LV Power Supply to operate the SiPMs

Both SiPMs are readout by means of a Lecroy HD06104 oscilloscope (1 GHz bandwidth, 2.5 GS/s*)

NB:

- SiPM 1 and SiPM 2 are operated at same gain (from previous measurements, $5x10^5 @ 56V$, ~room T)

- The trigger is active when both SiPM 1 and SiPM2 produce a signal ($V_{thr} > 5 \text{ mV}$) in a time window of 200 ns



Preliminary test: measuring the rise time



- We implemented an online analysis on the HDO6104 and collected the rise time of each signal.
- Working at very low thresholds exposes the system to EM noise.
- A filter might be used: the SNR will improve, however this would affect the shape of the signal, especially the fastest component.
- We decided to perform a more refined analysis at the waveform level.

Waveform analysis



The analysis is performed at the waveform level by measuring relevant parameters (Rise Time (20%–80%), Max Amplitude, Baseline)



Time delay: a good estimator for the time resolution



Considerations and conclusions

- The timing measurement of "scintillator bar + SiPM" is dominated by the instrumentation performance and the experimental setup. In this preliminary test we have collected SiPM signals in correspondence of crossing muons through a HD06104 Lecroy oscilloscope.
- The oscilloscope bandwidth (1 GHz) and sampling rate (2.5 GS/s) seem not to be adequate to perform a precise measurement.
- The absence of a preamplifier demands for a configuration in which the minimum discrimination threshold imposed is > 4 mV over 50 Ω @ a SiPM gain of 5 x 10⁵.
- 5 mV threshold corresponds to \sim 60 photons (and, due to the SiPM structure, the larger is the signal, the slower is the rise time).
- We could set an upper limit for the timing to about 400 ps (compatible with the actual sampling rate of the oscilloscope, 1 sample every 400 ps).

Room for improvement

- Move from analog to digital using Constant Fraction Discriminator and fats TDC
- Use identical scintillator bars, equalise finely the gain of the SiPMs, improve the optical coupling, clean the setup.
- Implement an external trigger logic to lower the threshold on the SiPM.