



FRIDA, INFN Pisa

Update on plastic scintillator dosimeter development (WP3)

E. Ciarrocchi (presenting), A. Cavalieri (Master student), D. Del Sarto, F. Di Martino, A. Kraan, M. Montefiori,

M. Morrocchi, J. Pensavalle, E. Ravera (PhD student), V. Rosso, M. Bisogni






Second GM, 29-30 March 2023

Quick recap

Plastic scintillator - Kuraray SCSF-78J, 10 mm length and 1 mm diameter

Clear optical fiber - Thorlabs, 1 mm diameter and 25 m long



Optical coupling - Cyanoacrylate glue in **Peek tube**






Performance of the first dosimeter. a) Linearity of the signal with the dose per pulse (DPP), for the two tested collimators, 100-mm diameter (blue) and 40-mm diameter (red). b) Dependence on the pulse repetition frequency (prf). c) Dependence on the intra-pulse dose rate (IPDR). d) Ratio of Cerenkov to total scintillator signal at different dose per pulse values.

Updates

- **Plastic scintillators (Kuraray)**
 - 10 mm length, 1 mm diameter (blue, 450 nm)
 - 2 mm length, 0.5 mm diameter (green, 530 nm)
- **Clear optical fiber** - 25/10 m long, matching the scintillator diameter
- **Optical coupling** - Thorlabs F120 optical glue in **carbon tube**
- Readout by a **CCD camera** (image of the fiber output) or by a fast **photodiode** (waveform)





PMMA/PVC block



1. Linearity of the blue scintillator



No monitor unit correction (see later)



General model:

$$\text{fitresult}(x) = a \cdot x / (1 + b \cdot x)$$



Coefficients (with 95% confidence bounds):

$$a = 0.057 \pm 0.006$$

$$b = 0.013 \pm 0.007$$


2. Percentage-depth-dose (PDD) curve

Reconstruction of the dose profile at variable depth in solid water




3. Beam spatial profile



Lateral beam profile at fixed depth in solid water



50-mm applicator




4. Beam time profile





5. Dosimeter aging

- We take a reference measurement @12 Gy/p
- We deliver several kGy
- We repeat
- Monitor unit correction applied!



+1 kGy approximately, due to previous irradiations



Future plans



- **Aging** of the plastic scintillator dosimeters @ CPFR and annealing with our X ray tube, by monitoring the scintillator output over time (minutes? days? weeks?)
- **Beam monitoring** with plastic scintillator sheet imaged by CCD camera
- Use of plastic scintillator dosimeters to measure the **entrance dose *in vivo*** with our colleagues of CNR Pisa

Papers and conferences in 2023

- Proffered paper at ESTRO 2023, Vienna, May 15th “Plastic scintillator-based dosimeters for FLASH radiotherapy”, session on Flash and proton measurements
- Paper in preparation for the Frontiers Special Issue “Multidisciplinary Approaches to The FLASH Radiotherapy”, deadline @end of April 2023