

Artificial Intelligence in Medicine



Working with clinicians: the DORIAN Tech. experience

Ruben Gianeri

next_AIM annual meeting

15.02.2023, Milano



DORIAN
evolving neuroimaging



DORIAN delivers a fast and reliable tool for the **quantification** of medical imaging to support the early and differential diagnosis of neurodegenerative disorders. It provides clinicians and researchers with state-of-the-art robust, rater-independent and reproducible quantitative biomarkers to **better evaluate** dementias stage and progression, **complementing** their ability to write informed medical reports and improve on the early detection and diagnosis of neurodegenerative diseases.

Dorian recent activity

- > Spin-off foundation [May 2020]
- > **Start of beta-tests @ 6 italian centers** [Apr 2021]
- > **Lab sessions @ Nucl. Med. school in Neurology** [Sep. 2021]
- > EBAN finalist [Mar.2022]
- > Unicredit StartLab 3rd classified [May 2022]
- >> Amazon AWS academy invitation [Sep. 2022]
- > **Lab sessions @ Nucl. Med. school in Neurology** [Oct. 2022]
- > First fulfilled contract [Dec. 2022]


A. Chincarini
(INFN-GE)



P. Bosco
(IRCCS Stella Maris, Pisa)



M. Corosu
(INFN-GE)




E. Peira
(INFN-GE)





R. Gianeri
(INFN-GE)



F. Sensi
(IRCCS San Martino, Genova)



Benefits of quantification



Inter-rater agreement of 6 expert clinicians in the differential assessment of 100 FDG-PET cases, which have been confirmed with a diagnosis of MCI-LB or MCI-AD.

The accuracy and agreement (w Cohen's k) of each clinician is plotted with respect to 3 different diagnostic settings:

1. "visual": PET data is presented as acquired,
2. "vis+spm" where the PET data is presented together with a consolidated semi-quantification algorithm
3. "vis+spm+quant" where the PET data and semi-quantification are complemented with a third machine learning-based method, and a comprehensive analysis is presented to the clinician.

The added information and its synthesis not only **improves the accuracy** of the single clinician with respect to the true diagnosis, but it also delivers a much **higher intra-rater agreement**.

There are several solutions for image analysis, but they are usually provided as workstations:

- > Difficult access to the workstation: the console is not necessarily in a space accessible to all clinicians, and it can only be used by one person at a time.
- > Rely on local computational resources.
- > Require constant maintenance and upgrading is usually done with the purchase of a new workstation.
- > Analysis are often restricted to certain Scanners only.
- > Semi-automated analyses: direct intervention of the clinician is required.



D0lab is

- Accessible from any device
- Computes on cloud
- Supports multiple users
- Easy to extend
- Secure, GDPR compliant
- Fully automated

1. **include clinicians/NM experts** in the process from the beginning
2. use of **multi-algorithm** approach, integrating standard methodologies with AI
3. aim at **quantification** (i.e. a direct product of the exam, open to interpretation) rather than at the diagnostic label (the outcome of a more complex process often involving several experts, difficult to challenge)
4. use **explainable** algorithms
5. **train users** & get feedback
6. keep implementation **simple**


Past activities

- > Collaboration with clinicians for early research development
- > Clinical testing @ 5 Italian centers: Genoa, Milan, Prato, Bergamo, Padua [2021 - 2022]
- > Training session at the AIMN National School of Nuclear Medicine (Advanced course) [Sep 2021]
- > Training session at the AIMN National School of Nuclear Medicine (Base course) [Oct 2022]

Planned activities

- > First release installation in AIMN-Neurology centers (15 Hospitals Nationwide) [May 2023]
- > First certified training course for DOlab usage in clinical practice [Dec 2023]

AIMN training [2021 & 2022]



Survey results

- about 50% of participants regularly use quantification tools.
- 92% would use DOLab in clinical practice.

Preferred features


- > multi-algorithm quantification
- > explainability of analysis algorithms
- > regional analysis
- > comparison of patient results with normative population

Absolute Scores Absolute Vs. Normative Scores Analysis Checks

ROI	SUVR		ELBA		TDr		CmpLS	
	Score	Cut-off	Score	Cut-off	Score	Cut-off	Order	Visual
L Frontal	1.85	1.15	1.09	0.8	0.68	0.54	1 st	●
L Occipital	1.23	1.2	0.68	0.71	0.49	0.47	9 th	●
L Posterior Parietal	1.47	1.24	0.9	0.81	0.64	0.55	3 rd	●
L Lateral Temporal	1.25	1.18	0.73	0.73	0.58	0.53	7 th	●
L Precuneus + Post Cingulate	1.56	1.32	0.86	0.9	0.58	0.58	8 th	●
R Frontal	1.8	1.23	1.08	0.81	0.65	0.57	2 nd	●
R Occipital	1.16	1.12	0.65	0.68	0.48	0.46	10 th	●
R Posterior Parietal	1.43	1.29	0.9	0.78	0.66	0.52	5 th	●
R Lateral Temporal	1.27	1.14	0.76	0.76	0.6	0.49	6 th	●
R Precuneus + Post Cingulate	1.55	1.25	0.88	0.83	0.6	0.53	4 th	●

CmpLS is around (●), below (●) or above (●) its cut-off

From research to clinical practice



Real world data

Heterogeneity of images:

- Quality
- Format
- Acquisition protocol



- Model validation on large multi-center dataset (~1000 clinically validated cases).
- Introduction of automatic checks on input data, providing warnings to the user.



Human mistakes

> Wrong image uploaded:

- NM image (i.e. Early AmyPET/FDG-PET instead of Late AmyPET)
- Wrong patient
- Flipped/preprocessed image

> Incomplete data

- Provide a comprehensive manual and adequate user training.
- Continuous checks and requests for confirmation of data submitted by the user.
- Automatic inspection on input data (FOV, resolution, ...).
- Image recognition with AI algorithms



Clinical Advisors



Diego Cecchin

Director of Nuclear medicine unit
@ University Hospital of Padua

“Many reliable analyses. A sound tool to increase diagnostic confidence!”



Silvia Morbelli

Nuclear Medicine Physician
@ IRCCS S.Martino, Genoa

“DOLab is a real gem! I use it everyday and I gladly suggest it to my colleagues!”




Stelvio Sestini

Director of Nuclear medicine unit
@ Prato Hospital

“Quantification made with the DOLab platform is really helpful to detect damaged brain regions and precisely monitor the therapy efficacy.”

Artificial Intelligence in Medicine



Thank you

Ruben Gianeri
next_AIM annual meeting
15.02.2023, Milano



DORIAN
evolving neuroimaging

