# HiDRa Simulation & Analysis Updates

15/02/2023

Andrea Pareti - INFN and Universita' di Pavia

Recap

HiDRa Sim with 80 mini-modules and 2mm outer capillary diameter

Pion resolution in [10, 100] GeV Range



Impact of Calo Geometry, from 84 to 120 and 480 mini-modules (Steel, 2mm capillaries)

Pion resolution in [10, 100] GeV Range



#### Recap

Fixing Steel absorber and 480 miniM, change capillary outer radius

#### Pion resolution in [10, 100] GeV Range



#### Results from Marco Lucchini (INFN & UniMiB)



Recap

Fixing Steel absorber and 480 miniM, change capillary outer radius

Pion resolution in [10, 100] GeV Range



#### Results from Marco Lucchini (INFN & UniMiB)





• Implemented Marco's Fiber disposition



- Implemented Marco's Fiber disposition
- Studied Optical photons propagation inside the Cerenkov fibers.

Marco's Sim doesn't consider

- Possonian Smearing Effect
- Total Reflection inside fibers



#### Cerenkov Phe/GeV dependence on energy

#### Standard HiDRa, without Smearing and Total Reflection

Pion resolution in [10, 100] GeV Range, Linear Error sum



Our "Best case scenario": brass absorber, 480miniM, alternative fiber disposition and nor Smearing neither Total Reflection

Pion resolution in [10, 100] GeV Range, Linear Error sum



### • Implemented Marco's Fiber disposition

 Studied Optical photons propagation inside the Cerenkov fibers.

Marco's Sim doesn't consider

- Possonian Smearing Effect
- Total Reflection inside fibers
- ROOT Fit extracts different sampling and constant terms, but resolutions are now closer

#### **Pion Resolution Comparison**



### • Implemented Marco's Fiber disposition

 Studied Optical photons propagation inside the Cerenkov fibers.

Marco's Sim doesn't consider

- Possonian Smearing Effect
- Total Reflection inside fibers
- ROOT Fit extracts different sampling and constant terms, but resolutions are now closer

Useful exercise to understand the importance of optical effects in the simulation, and of the Cerenkov photon yield

#### **Pion Resolution Comparison**

