



UPDATE ON NIT EMULSION PILOT TEST



FragementatiOn
Of Target

A. Alexandrov, T. Asada, V. Boccia, N. D'Ambrosio, A. Di Crescenzo, G. De Lellis, G. Galati, A. Iuliano, A. Lauria, G. Lorusso, M.C. Montesi, V. Tioukov

Università di Napoli "Federico II", INFN Napoli
Università di Bari "Aldo Moro", INFN Bari
INFN LNGS

05/06/2023, XIV FOOT General Meeting, Bergamo

PRIN Proposal Approved!

Project title: **DAMON: Direct meASureMent of target fragmentatiON**

Coordinator: **GALATI Giuliana** ERC: **PE2_6**



University: **Università degli Studi di BARI ALDO MORO**

n°	Associated Investigator	Qualification	University/ Research Institution
1.	GALATI Giuliana	Ricercatore a t.d. - t.pieno (art. 24 c.3-a L. 240/10)	Università degli Studi di BARI ALDO MORO
2.	LAURIA Adele	Professore Associato (L. 240/10)	Università degli Studi di Napoli Federico II
3.	D'AMBROSIO Nicola	Dirigente tecnologo	Istituto Nazionale di Fisica Nucleare





Direct Measurement of Target Fragmentation

- The goal of the project is the *direct* measurement of target fragments produced by a proton beam
- Nano Imaging Trackers (NIT) emulsions act both as target and tracking devices
- Each NIT film has two sensitive layers (40 μm thick) deposited on both sides of a plastic support (200 μm thick)






Data Taking in Trento

- 19 NIT films ($\sim 6.4 \times 4 \text{ cm}^2$) for the brick, 1 film for sensitivity tests
- Fixed pencil beam @211 MeV (FWHM $\sim 1.5 \text{ mm}$)
- 230 cm from beam exit window
- 6x4 grid for a uniform exposure of the NIT emulsions (about 11.000 protons per spot)




Data Taking in Trento

- The distance from the beam exit window was optimized to achieve a uniform exposure of the brick with TOPAS simulations
- Goal: 10.000 protons/cm²
- First spot with about 85.500 protons instead of 11.000
- Two spots (20.000 and 40.000 protons) for sensitivity tests





First Scanning Tests

- A preliminary scan was performed in LNGS on the top side of R1-2
- 100x100 views, each view being $50\ \mu\text{m} \times 70\ \mu\text{m}$
- Assuming a fluence of $10.000\ \text{protons}/\text{cm}^2$ about 3500 primaries are expected (about 1 in every 30 views)





Actual emulsion layer thickness around $35\ \mu\text{m}$



Micro-track Reconstruction


- Once grains are reconstructed, the linking procedure takes place
- The scan performed @LNGS was used to optimize the parameters for the linking algorithm
- Among the parameters:
 - Minimum Number of grains = 2
 - Maximum kink angle between grains = 0.8 radians
 - Maximum linking distance = 5 μm




Micro-track Reconstruction

- A score can be given to each microtrack: $Score = -\log(\langle gap \rangle \cdot \langle kink \rangle)$
- Cuts on the score and the number of grains identify the best microtracks
 - At least 10 grains, score > 0.75 ($\langle gap \rangle \sim 1 \mu m$ and $\langle kink \rangle \sim 30^\circ$)
- The horizontal micro-tracks are closer to the surfaces and they in part due to scratches
- The vertical micro-tracks are less than expected


Average Distance between micro-track grains



Micro-tracks Angular Distribution



Micro-track θ vs Z




Microscope Upgrades

- Target fragmentation measurements with NIT emulsions → need new hardware!


Standard FOOT microscope

- Standard nuclear emulsions (grain size = 200 nm)
- Works in transmission
- 20x objective
- 7 μm physical pixel size → 350 nm pixel size
- Z Step = 1.75 μm
- High scanning speed (~ 20 cm^2/h , up to 190 cm^2/h)
- Small grains **not** visible



- NIT emulsions (grain size down to 40 nm)
- Works in reflection
- 100x objective + 2.6 zoom lens
- About 27 nm pixel size
- Z Step = 250 nm
- Slow scanning speed (about 14 days to scan top layer of one exposed film! ~ 0.07 cm^2/h)
- Small working distance, only one side at a time
- Blue light → plasmonic resonance

NEWSdm microscope





New Microscope to scan NIT emulsions

- Our needs are intermediate between the existing solutions → new design
- New optics arrived in May
- Microscope assembly almost completed!

New FOOT microscope



- Works in reflection
- 20x/40x objective + **adjustable** magnifying lens
- Intermediate Z step
- Better scanning speed
- Larger working distance, possible to scan both sides
- Blue light → plasmonic resonance



New Microscope Tests (1)

- First images with 20x objective (NA=0.75) installed (no magnification lens)
- Light source: $3 \times 4 \text{ mm}^2$ blue LED ($\lambda = 450 - 470 \text{ nm}$)

Field of View: $\sim 800 \times 600 \mu\text{m}^2$



New Microscope Tests (2)


- First images with 20x objective (NA=0.75) installed (no magnification lens)
- NIT sample exposed to Am alpha source (Trento exposure batch)
- Good contrast for alpha tracks




New Microscope Tests (3)

- Comparison between 20x (NA=0.75) and 40x objective (NA=1.3)
- R1-17 film from the brick
- Manual search for interaction vertices

$$\text{Rayleigh criterium: } R = 0.61 \frac{\lambda}{NA}$$



40x objective



20x objective

Conclusions

- Preliminary scan @LNGS → optimization of linking algorithm parameters
- Assembly of new microscope almost completed (ETA ~ 1 month)
 - Correct light source position
 - Optimize Z step
 - Test adjustable magnifying lens
 - Test top side scanning / bottom side scanning
- Good contrast observed for alpha tracks (NIT sample poured on glass)
- Difficulty in reconstructing vertical proton tracks!
 - To be confirmed with full area scans: further studies on-going