

TIME-RESOLVED PHOTOEMISSION EXPERIMENTS WITH LASER BASED HHG SOURCES

Barbara Ressel

Laboratory of Quantum Optics
University of Nova Gorica, Slovenia

Frascati, 6.12.2024



OUTLINE

CITIUS light source

High order Harmonic Generation (HHG)


Applications:

Time resolved photoemission


Challenges and examples



CITIUS LIGHT SOURCE



HHG GENERATION



Possibility of generating lower harmonics using nonlinear crystals:

- 2nd (400nm) PE = 3 eV
- 3rd (266nm) PE = 4.5 eV
- 4th (200 nm) PE = 6 eV

FERMI FEL-1


Energy range	12.4 — 120 eV
Pulse duration	20 — 40 fs
Rep rate	10/50 Hz
Photons/pulse	$10^{12} - 10^{14}$

CITIUS




Energy range	14 — 80 eV
Pulse duration	35 fs
Rep rate	5 kHz
Photons/pulse	$10^5 - 10^8$



TIME RESOLVED PHOTOEMISSION



TYPICAL DATA SET




CHALLENGES

- ✓ LIMITED ENERGY RESOLUTION
- ✓ REPETITION RATE
- ✓ MULTIPHOTON ELECTRON EMISSION
- ✓ SPACE CHARGE EFFECT
- ✓ LASER ASSISTED PHOTOEMISSION (LAPE)



MULTIPHOTON EMISSION




CAUSES

-THERMIONIC EMISSION or HOT EMISSION: just like emission from a hot filament, but the electrons are heated with the laser.



- FIELD EMISSION or COLD EMISSION: the potential of the atom is perturbed by the laser's electric field so that the electrons are ripped from the parent atom.






EFFECTS

- BACKGROUND
- ENERGY SPREAD AND SHIFTS (due to the resulting space charge effect)


REMEDIES

- Reduction of laser energy flux
- Alignment and sample preparation to avoid illuminating facets and defects

SPACE CHARGE EFFECT






LASER ASSISTED PHOTOEMISSION



In a LAPE process additional photon can be absorbed or emitted by the electron, resulting in sidebands.


LAPE can be very useful as a tool for the characterization of XUV pulse duration.



Can form replicas of spectral features in a PES experiment.

Can induce features not related to the response of the material.

ENERGY RESOLUTION



PtTe₂ ARPES data - Courtesy of dr. Zacchigna (CNR-IOM)



THE “GOOD” SIDE

- ✓ Straight forward control on the linear polarization
- ✓ Possibility of further compression of pulses (hollow fiber compressor)
- ✓ GAS PHASE PHOTOEMISSION - ref. dr. Marcello CORENO



SELECTIVE ORBITAL EXCITATION

Both pump (IR) and probe (XUV) beams allow for easy manipulation of their linear polarizations



PHYSICAL REVIEW B **98**, 205142 (2018)

Orbital-dependent electron dynamics in Fe-pnictide superconductors




Ganesh Adhikary,^{1,*} Barbara Ressel,¹ Matija Stupar,¹ Primož Rebernik Ribič,² Jurij Urbančič,¹ Giovanni De Ninno,^{1,2,†} D. Krizmancic,³ A. Thamizhavel,⁴ and Kalobaran Maiti^{4,‡}

Using polarization to selectively excite orbitals in EuFe_2As_2 pnictide

SELECTIVE ORBITAL EXCITATION



Europium pnictide - EuFe_2As_2

- ✓ Parent compound to Fe superconductors
- ✓ Study of electron correlation







Different dynamics for different symmetry bands are fingerprint of different interactions and localizations of bands

HOLLOW FIBER COMPRESSOR



ULTRAFAST OSCILLATIONS IN TaS₂



Simoncig A., et al., Phys. Rev. B 103, 155120

CONCLUSIONS

CITIUS novel light source providing ultrashort light pulses in the XUV soft X ray spectral region

TR photoemission pros and cons

Access through proposal submission

- NEP (NFFA Europe Pilot), call every 3 months

<https://www.nffa.eu>



- Elettra Sincrotrone Trieste - Fermi, call every 6 months

<https://www.elettra.eu/userarea/access-request.html>



Elettra Sincrotrone Trieste



LKO GROUP AND COLLABORATORS

LKO @ University of Nova Gorica

Prof. Dr. Giovanni De Ninno

Izr. Prof. Dr. Barbara Ressel

Dr. Johannes Schmidt

Ph.D. students:

Rokaya Osama

Arun Ravindran

Luka Novinec

Sincrotrone Trieste S.C.p.A.

Prof. Dr. Giovanni De Ninno

Dr. Primož Ribič Rebernik

Dr. Carlo Spezzani

Luka Novinec

Dr. Alberto Simoncig

Dr. Marcello Coreno (ISM-CNR)

CNR- IFN (Padova, Italy)

Dr. Luca Poletto

Dr. Fabio Frassetto

Dr. Paolo Miotti

