

# Progress of an accelerator mass spectrometry system on the Tsukuba 12UD Pelletron tandem accelerator

-11th International Conference on Heavy Ion Accelerator Technology - Venezia (Italy) –  
8-12, June, 2009

Kimikazu SASA

Tandem Accelerator Complex (UTTAC),  
University of Tsukuba, Japan

## Collaborators


T. Takahashi, Y. Nagashima, Y. Tosaki, K. Sueki,  
Y. Yamato, N. Kinoshita, T. Amano,  
UTTAC, University of Tsukuba  
H. Matsumura, K. Bessho, RSC, KEK  
Y. Matsushi, MALT, The University of Tokyo



# Outline of presentation

---

- Introduction
  - 12UD Pelletron tandem at the University of Tsukuba
  - AMS and facilities
- AMS system on the 12UD Pelletron tandem
  - Description of the Tsukuba AMS system
  - Recent progress
  - Performance of  $^{26}\text{Al}$ ,  $^{36}\text{Cl}$  and  $^{129}\text{I}$  AMS
- Summary and future plans



University of Tsukuba,  
Tsukuba science city

60 km from Tokyo in the northeast

2 accelerator facilities


Proton Medical Research Center: PMRC

Tandem Accelerator Complex: UTTAC



• 250 MeV Proton Synchrotron (2001)

➔ Proton Beam Radiotherapy



• 12UD Pelletron Tandem Accelerator (1975)  
• 1MV Tandatron Accelerator (1987)



# 12UD Pelletron tandem accelerator (1975)

University of Tsukuba, Tsukuba, Japan.


- MODEL : Vertical Tandem Van de Graaff
- Terminal voltage range: 1 - 12 MV
- Insulation Gas : SF<sub>6</sub> pressure: 0.6 MPa
- Accelerator Tank : Height: 17.9 m  
Diameter: 4.8 m  
Volume: 350 m<sup>3</sup>  
Total weight: 120 ton
- Analyzing Magnet : ME/Z<sup>2</sup>=200 amu MeV

In the early stage


Nuclear physics

Now

- AMS
- IBA
- Material science



# 12UD Pelletron tandem accelerator (1975)



# Upgrade of the 12UD Pelletron tandem



## 2009 Divided resistor system

We replaced the old corona needles with the divided resistor system.


Variable terminal voltage  
(No shorting column) →

$$V_t = 1 \sim 12 \text{ MV}$$

# Accelerator Mass Spectrometry

## Targets of AMS

$^{10}\text{Be}$  ( $T_{1/2}=1.36\times 10^6$  yr)  
 $^{14}\text{C}$  (5730 yr)  
 $^{26}\text{Al}$  ( $7.1\times 10^5$  yr)  
 $^{36}\text{Cl}$  ( $3.0\times 10^5$  yr)  
 $^{129}\text{I}$  ( $1.57\times 10^7$  yr)  
...




Isobar, Isobaric Molecular suppression

Accelerator



Mass analysis



Detectors

High energy part


Momentum analysis

energy analysis

$\sim 1$  MeV/amu

# Beam energy for the 12UD Pelletron

Terminal voltage: 1 - 12 MV



- AMS on the 12UD Pelletron tandem
  - Description of the Tsukuba AMS system
  - Recent progress
  - Performance of  $^{26}\text{Al}$ ,  $^{36}\text{Cl}$  and  $^{129}\text{I}$  AMS

# Tsukuba AMS system

## Progress of the Tsukuba AMS system

1993-1996	Trial AMS measurement for $^{14}\text{C}$ .
1996-1998	Development of AMS system $^{14}\text{C}$ -AMS AMS ion source (original) Mass separator beam line
1999-	Development of $^{26}\text{Al}$ , $^{36}\text{Cl}$ -AMS Pilot beam methods (Instead of GVM control)
2002-	Development of $^{129}\text{I}$ -AMS
2007-	$^{36}\text{Cl}$ AMS 9 MV $\rightarrow$ 10 MV (Improved beam line) Background: $^{36}\text{Cl}/\text{Cl} < 1 \times 10^{-15}$ Repetition accuracy: $\pm 3\%$
2009-	Upgrade of the 12UD Pelletron (Resister system)  GVM terminal control system

# Tsukuba AMS system

A pilot beam method is used to stabilize the terminal voltage.


➔ The terminal voltage is kept stable within 0.02 %.

$^{26}\text{Al}$ ,  $^{36}\text{Cl}$ ,  $^{129}\text{I}$  - AMS


Pilot beams (isobar)



500 samples/year.




# $^{36}\text{Cl}$ -AMS by the Tsukuba AMS system




Target material	AgCl + C <sub>60</sub>
Cl <sup>-</sup>	~ 20 mA
V <sub>T</sub>	10 MV
Pilot beam	$^{12}\text{C}_3^-$
Detection ion	$^{36}\text{Cl}^{14+}$
Particle energy	100 MeV
Detection range	$^{36}\text{Cl}/\text{Cl} = 10^{-10} \sim 10^{-14}$
Back-ground	$^{36}\text{Cl}/\text{Cl} < 10^{-15}$
Repetition accuracy	± 3 % ( $^{36}\text{Cl}/\text{Cl} \sim 10^{-12}$ )

# AMS Cs sputtering ion source




# Upgrade of LEBT for 12UD Pelletron tandem accelerator



Einzel lens was replaced with electrostatic quadrupole (triplet).

Beam transport  
Injection energy  
 $^{36}\text{Cl}^-$  : 103 kV  
H $^-$ : 120 kV  
Terminal voltage: 10 MV




# Tsukuba 12UD first unit (2009)

Lower resistance  
at the entrance.




Adjustment for the focusing effect




# Terminal section (Charge exchange)

Terminal section was modified to the large aperture canal ( $\phi 20$ ) in 2004.




$\phi 16$   
for AMS


Carbon foil for AMS:  $5 \text{ mg/cm}^2$



# Mass separator beam line




# Gas E- $\Delta E$ detector



100 MeV  $^{36}\text{Cl}$  in the gas detector.

# $^{36}\text{Cl}$ -AMS (2-dimensional spectrum)




Standard sample

$^{36}\text{Cl}/\text{Cl} = 1.60 \times 10^{-12}$

Blank sample

# $^{26}\text{Al}$ -AMS (2-dimensional spectrum)



- Full stripping technique

- Pilot beam:  $^{26}\text{MgO}^-$

- Beam current of  $\text{AlO}^-$  from  $\text{Al}_2\text{O}_3$  sample :  $\sim 1.5$  mA

-  $^{26}\text{Al}$  is very clearly separated from  $^{26}\text{Mg}$ .

- Background of the  $^{26}\text{Al}$ -AMS:  $< 1 \times 10^{-15}$ .

# Performance of the Tsukuba AMS system

A pilot beam method is used to stabilize the terminal voltage.

## $^{26}\text{Al}$ -AMS

Target material	$V_T$	Injection ion	Pilot beam	Detection ion	Particle energy	Back-ground
$\text{Al}_2\text{O}_3 + ^{26}\text{MgO}_2 + \text{Ag}$	10.2 MV	$^{26}\text{AlO}^-$	$^{26}\text{MgO}^-$	$^{26}\text{Al}^{13+}$	78 MeV	$< 1 \times 10^{-15}$


## $^{36}\text{Cl}$ -AMS

Target material	$V_T$	Injection ion	Pilot beam	Detection ion	Particle energy	Back-ground
$\text{AgCl} + \text{C}_{60}$	10 MV	$^{36}\text{Cl}^-$	$^{12}\text{C}_3^-$	$^{36}\text{Cl}^{14+}$	100 MeV	$< 1 \times 10^{-15}$


## $^{129}\text{I}$ -AMS

Target material	$V_T$	Injection ion	Pilot beam	Detection ion	Particle energy	Back-ground
$\text{AgI} + \text{MoO}_2 + \text{Nb}$	9.7 MV	$^{129}\text{I}^-$	$^{97}\text{MoO}_2^-$	$^{129}\text{I}^{26+}$	126 MeV	$< 1 \times 10^{-13}$

# Applications by the Tsukuba AMS system


 Mainly for earth and environmental sciences.

**Nuclear safety research**  
 Atomic bomb, neutron fluence



Hiroshima A-bomb sample

**Soil sediment**



Soil

**Rock meteorite**



Limestone

Meteorite

**Groundwater, rain, ice**



Rain water



Ice core

**Biological sample**



Human hair

# Summary and future plans

- 12UD Pelletron tandem at the University of Tsukuba

We have upgraded the 12UD Pelletron tandem.

⇒ LEBT, divided resistor system, terminal stripper.

The beam time for AMS research has increased to about 42% of the total operation time.

- Tsukuba AMS system

We are able to measure long-lived radioisotopes of  $^{26}\text{Al}$ ,  $^{36}\text{Cl}$  and  $^{129}\text{I}$  by employing a molecular pilot beam method that stabilize the terminal voltage with 0.02% accuracy.

⇒ Main research fields are earth and environmental sciences.

## Future plans

- GVM control system

- New injection beam line (MC-SNICS)

Thank you for your kind attention.

