CHNet_ILAGING Development of a transmission imaging system for muon applications*

Matteo Cataldo, *PhD Proposta per call «Grant Giovani 2025»*29/07/2025

*acronimo e titolo da definire

About me: research fellow @ INFN Milano Bicocca

My Career:

- 2014 2019: Science and Conservation for Cultural Heritage
 @Università degli studi di Firenze
- 2020 2024: PhD @Università degli studi di Milano Bicocca/<u>ISIS</u> <u>Neutron and Muon Source</u>
- 2024 Present: Research fellow @INFN Milano Bicocca



ISIS Neutron and Muon Source





Development of the Muonic atom X-ray Emission Spectoscopy (μ-XES) technique

+


Use of large scale facilities techniques for the characterisation of Heritage artefacts

+

CHNet_MAXI (FTE 0.6)

-


RESNOVA (FTE 0.4)



With my PhD supervisors, Massimiliano Clemenza and Adrian D. Hillier

Introduction: negative muon spectroscopy

Muonic atom X-ray Emission spectroscopy (μ-XES) is the technique based on the detection of high energy radiation (X and gamma ray) emitted after muon interaction with matter




Multi-elemental

Non-invasive probe
No activation of samples

Method for isotopic analysis

Negligible selfabsorption effect

Remarkable penetration depth Momentum & Density dependent











Introduction: bringing cosmic-ray muon imaging to the laboratory scale


Cosmic muons are used in two ways^[1]:

Muography (Muon trasmission radiography)

- Detector is placed below or next to the object
- The attenuation of the muon flux gives information about the density of the material
- Observed muon flux vs open-sky flux to get a density map


Muon Scattering Tomography (MST)

- Muons undergo multiple <u>Coulomb scattering</u> when travelling through matter. **Scattering depends on:**
 - Atomic number
 - Thickness and density of the material.
 - Muon energy
- Collection of muon tracks (in&out with a tracking detector) to produce a **scattering map**



Is it possible to scale down these applications to a muon facility?

Introduction: bringing cosmic-ray muon imaging to the laboratory scale




Applications of muography and MST:

Archaeology, geosciences, civil engeneering, nuclear waste analysis, border control

We need to scale down the samples too!





Goal of the project

Develop a high-resolution, portable muon imaging setup to support Muonic Atom X-ray Spectroscopy (μ-XES) experiments and enable independent, non-destructive imaging.

Where?

INFN Milano Bicocca

+

ISIS Neutron and Muon Source

(Controlled muon beam <u>up to</u> 120 MeV/c)

How?

Experimental setup: Scintillator screen


+

CMOS Camera to record the light


Why?

- 1) Support and improve μ-XES capabilities, to make it a unique elemental analysis method;
- 2) Explore the possibilities offered by muon imaging
- 3) Continue INFN experience in developing large scale facility techniques (CHNet Tandem, Niche, Bronze, MAXI)

Literature review



@J-PARC: Dose range estimation in water and plastic scintillator



@J-PARC: Development of a beam profile monitor

Proposed experimental setup

Positive and/or negative


Incident muon

beam


Transmission measurement to produce a muon


radiography











Light

GAGG/Plastic/LYSO scintillator (to be evaluated)




Detector



Focusing lens + **CMOS** sensor

Can we get through?

Consider a 5mm thick silver coin (Ag90Cu10 – Density 10.3 g/cm3). A 80 MeV/c muon beam will travel the entire size of the coin.



Expected outcome

- Enhanced muon experiment accuracy through imaging
- Safe, non-destructive, activation-free analysis
- Proof-of-concept for broader muon imaging applications
- Tool for better beam targeting in complex samples or holders
- Estimated sub-mm resolution



Possible application fields:

- Heritage Science
- Battery materials
- Defects identification in materials
- Wood, Paper material?

Pros and Cons

Pros

- No activation of sample
- Good penetration depth
- Positive and negative muons can be used
- No radiation damage on sample
- Spatial resolution: sub-mm (in literature, similar setup 300-400 nm)
- High muon flux
- Multisciplinary
- Innovative

Cons

- <u>Competitors</u>: Neutron and X-ray tomography
- Maximum momentum 120 MeV/c
- Image blurring due Coulomb scattering.
- Scintillator may need to be close to the sample to have a good rate.
- Muon beam spot size (nominal diameter 40mm)
- Few literature examples

Working packages

- WP1: Scintillators characterisations. Light output, linearity, decay time, and afterglow test to be performed at the Milano Bicocca. Test will be performed with different scintillator thickness (as it is a trivial parameter for resolution). Modelling of the setup in Geant4. Building of the light-tight box and the XY stage in collaboration with the INFN workshop. The light-tight box will be 3D printed.
- WP2: CMOS Camera characterisation. In house test of the setup with known sources and the CMOS camera in place. Evaluation of camera properties. Optimization of setup: mirror and without mirror test. Pinhole configuration test.
- WP2: First beam acquisition. Evaluation of the characteristics of the scintillator under a muon beam.
 Radiation damage evaluation on the scintillator foils. Test with positive and negative muons to evaluate the byproduct of the interaction. This test will be performed at the ISIS Neutron and Muon source. Preliminary data treatment
- **WP4: Final beam test.** Analysis of standard materials. Development of a data analysis process (MLEM/POCA algorithms?) and implementation/test of a rotating stage-

WP1 - Scintillator material

Inorganic	GAGG(Ce)	LYSO(Ce)
Density (g/cm3)	6.6	7.25
Light output (ph/MeV)	From 30000 to 54000 (depending on type)	30000
Decay constant (ns)	From 50 to 150 (depending on type)	40
Emission wavelenght (nm)	<u>520-530</u>	420
Hardness (mohs)	8	5.9
Hygroscopic	No	No

Plastic	(Epic 106)	EJ-212
Density (g/cm3)	1.05	1.023
Light output (% Anthracene)	50-60	65
Light output (ph/MeV)	-	10000
Decay constant (ns)	2.4	2.4
Emission wavelenght (nm)	423	423
Hygroscopic	No	No


GAGG: high light yield, emisssion wavelenght in the best working range of CMOS, radiation hard. But more expensive (1mm 8x8:1.4k€)

LYSO: similar to GAGG, but less exspensive: (1mm 8x8:900€).

Plastic: low density so thickness is an important choice. Quite cheap (1mm 8x8:150€)

Pixelated array or screen?

WP1 - Scintillator: pixelated or screen?



Screen: Image blurring due to light scattering, potentially reducing spatial resolution.

Low cost and thin size

Pixelated: reduce light scattering for better spatial resolution.

Is it possible to have low thickness?

High cost?

WP2 – CMOS Camera

ZWO ASI2600MC MONO

Pros:

- Extremely cheap (2k€)
- Easy to use
- Easy to adapt with focal lenses
- Already used in Niche
- Good resolution
- High pixel density
- C-mount
- Pixel 3.75 μm

Cons:

- No possibility of external triggering
- Low frame rate (3.5 fps)



Hamamatsu Orca Fusion

on

Pros:

- Good low-light performance
- Low read-noise
- High frame rate (up to 89 fps)
- External trigger
- Good resolution
- Pixel binning is feasible
- C-mount
- Pixel 6.5 μm

Cons:

- High cost (about 17 k€)
- Requires vendor software



Milestones and people

- Milestone 1 (month 6): From a set of different types of scintillators, choice of the most appropriate one in terms of size, thickness and light yield.
- Milestone 2 (month12): Setup optimization. CMOS Camera test and preliminary approach to data processing.
- Milestone 3 (month 18): First beam test. Evaluation of the scintillator performance under a muon beam and evaluation of radiation damage and by-product of the muon interaction.
- Milestone 4 (month 24): Final test of standard samples and a case study.

People	Sezione	WP	Fte	
Matteo Cataldo	Mib	1 to 4	0.8	
Paola Monza	Mib	1-2	0.4	
Giulia Marcucci	Mib	3-4	0.4	
Massimo	Pv	1	0.1	
Rossella				
Massimiliano	Mib	1 to 4	0.2	
Clemenza				
Luisa Vigorelli	a Vigorelli Mib		0.4	
			Total = 2.3	
Officina	Officina Mib		1MU	

External support	Institution
Adrian Hillier*	ISIS Neutron and Muon Source
Tom Ashton*	UCL
Matt Wilson*	ISIS Neutron and Muon Source
*Will write a support letter	

Funding request

	Detailed costs table	
	Plastic scintillator	1,500.00€
Consumables	GAGG	2,000.00€
	LYSO	1,500.00€
	Mirror	100.00€
	RS travel case	500.00 €
	Hammamatsu Fusion CMOS Camera*	17,000.00 €
Inventoriables	Acquisitoin Laptop	2,500.00€
	ASI2600MM-PRO CMOS Camera (Mono)*	2,000.00€
	Camera accessories	1,000.00€
	Experimental Activity** (for at least 4 runs for 2 people and eventual run in other facility)	10,000.00 €
Travel expsenses	Meeting	2,000.00€
	Conferences (at least one confernence)	2,000.00 €
	Total	<mark>42,100.00 €</mark>
still have two decide between the two. Or fpossible to buy both		
**Estimated cost for 3 days of beam time @RAL	3 nights in guesthouse+flights+forfeit+trains	880.00€

Still thinking on scintillator material, type and number of pieces. Also, we may need to ass some electronic equipment and cables

Estimated 50 k€ (?) Dividere in due anni

Risk Evaluation

Risk Chart		Impact▶					
	iok Oriar	_	Insignificant 1	Minor 2	Moderate 3	Major 4	Catastrophic 5
À	Almost Certain (>1 in 10)	5	Medium 5	High 10	Extreme 15	Catastrophic 20	Catastrophic 25
 pc	Likely (1 in 20)	4	Low 4	Medium 8	High 12	Extreme 16	Catastrophic 20
Likelihood	Occasional (1 in 200)	3	Low 3	Medium 6	High 9	High 12	Extreme 15
Lik	Unlikely (1 in 2000)	2	Very Low 2	Low 4	Medium 6	Medium 8	High 10
	Rare (<1 in 10000)	1	Very Low 1	Very Low 2	Low 3	Low 4	Medium 5

Item	Risk	LIKELIHOOD	Effect	IMPACT	RISK SCORE	Mitigation strategies
Scintillator	Delivery delay	3	Delay on WP1 start	3	9	The ISIS Facility will made availabe a set of spare scintillators if we experience major delays
Scintillator/ mirror	Radiation damages	4	Material need to be replaces	2	8	Test will be performed to evaluate this aspect. If present, it could be mitigated by fluxing argon in the light-tight box and by having spare scintillating material
CMOS Camera	Delivery delay	2	Delay on WP3 start	2	4	WP3 is expected to start in the firt half of 2026, so delay can be handled
ISIS Neutron and Muon Source	Proposal refusal	2	No beam test	5	10	As a group, we have a long experience of accepted proposal, so it's unlikely to be refused
ISIS Neutron and Muon Source	Shut down for entire 2027	4	No beam time in 2027, WP4 at risk	5	20	The facility refurbishment is expected to finish in August 2027, which leaves 2 run of beam time. To mitigate, we could submite proposals to other facilities
Material	Increase in cost	3	increase in cost could affect the funding request	3	9	a 10% will be added to the exsisitng quotes to cover increase in costs

Backup - Quotes

INFN MILANO BICOCCA MATTEO CATALDO P.ZZA DELLA SCIENZA, 3 20126 MILANO (MI) ITALY

Ns. Riferimento:

Tel:

Fax:

E-Mail:

Offerta n.: 1117478 11-07-2025 Data preventivo: Revisione n.:

HPI000651

Codice cliente: MATTEO CATALDO

Vs. Rif.: RFQ n.:

Pagina 1/1

francesco.giovenale@hamamatsu.eu

Francesco Giovenale

00390650513454

00390293581741

Pos.	Codice prodotto Descrizione	Quantità (pcs)	Prezzo unit. EUR	Importo netto EUR
1.1	ORCA-FLASH4 V3 USB FUL ORCA-Flash4 V3 with USB3 frame grabber and cable, trigger cable, HCImage Live ROHS:C	1	13.500,00	13.500,00
2.1	ORCA-FUSION USB FUL ORCA-Fusion with USB3 frame grabber and cable, trigger cable, HCImage live ROHS:C	1	14.300,00	14.300,00

RoHS: C = Conforme; L = Può essere utilizzato solo per dispositivi medici, strumenti di controllo e misura; N = Non conforme; U = Contatta l'ufficio vendite Hamamatsu locale

11-09-2025 Validità offerta: Porto Franco Termini di spedizione: Termini di pagamento: R.D. 30 GG. D.F. 12 Mesi Garanzia: 45 GG D.R.O. Tempi di consegna:

I prezzi si intendono al netto dell'aliquota IVA.

ARTESKY SRL

VIA IV NOVEMBRE 103, 20833 Giussano (MB)

Telefono 0362310657

AMMINISTRAZIONE@ARTESKY.IT Email

PEC ARTESKYSRL@PEC.IT

12059300967 P.IVA 12059300967

CCIAA e Numero REA: MB 2637730 Capitale Sociale i.v.: € 10.000,00

Spettabile

INFN ISTITUTO NAZIONALE DI ASTROFISICA SEZ MILANO

VIA CELORIA 16

20133 Milano (MI)

C.F. 84001850589 P.IVA 84001850589

VS EMAIL DEL 20/06/2025 Riferimento

DESCRIZIONE	QUAN	ITITA	IMPORTO	SCONTO	IVA	TOTALE
ASI2600MC-PRO	1,00	pz	1.589,00 €	7	22%	1.477,77 €
ASI2600MM-PRO	1,00	pz	2.120,00 €	7	22%	1.971,60 €
ASI6200MC-PRO	1,00	pz	3.950,00 €	10	22%	3.555,00 €
ASI6200MM-PRO	1,00	pz	3.950,00 €	7	22%	3.673,50 €
Termini di pagamento					Imponibile	10.677,87 €

rermini ai pagamento

Bonifico Modalità:

30 gg Data Fattura Scadenze:

Banca mittente

Banca: INTESA SANPAOLO SPA IBAN: IT97P0306951371100000010311

2.349,13 € Imposta 22% 13.027,00 € Totale

> 10.677,87 € Netto da pagare