Nuclear emulsions in FIRST: a first look




Giovanni De Lellis, Lucia Consiglio, Adele Lauria, Maria Cristina Montesi University Federico II and INFN, Naples, Italy

Nuclear emulsion as sensitive media for charged particles

After the passage of charged particles through the emulsion layer, a latent image is produced

After a physical/chemical process known as development, Ag grains become visible with an optical microscope



AgBr crystal, size 0.2 micron Is the elementary detection element

Nuclear emulsions used for more than 100 years in Particle Physics

Recorded as silver grains along the line particle passed through

30-40 grains/100 microns for MIP

50 micron

Resolution of 0.3 micron

01/02/2013


Microscopic Image

OPERA emulsions

OPERA industrial emulsions from FujiFilm

- The AgBr density in the OPERA emulsions is higher in respect to the commercial films
- •Special R&D for OPERA: the double pouring procedure

Emulsions are continuously sensitive detector ALL charged particle: cosmic rays, natural radioactivity etc recorded as a latent images. They can be partially erased by a "refreshing" procedure applied just before the detector assembling




Emulsion Scanning Stations

EU: ESS (European Scanning System)



- Scanning speed/system: 20cm²/h
- Customized commercial optics and mechanics
- Asynchronous DAQ software

Japan: SUTS (Super Ultra Track Selector)



- Scanning speed/system: 75cm²/h
- High speed CCD camera (3 kHz),
 Piezo-controlled objective lens
- FPGA Hard-coded algorithms

Both systems demonstrate:


- ~0.3 µm spatial resolution
- ~2 mrad angular resolution
- ~95% base track detection efficiency

Z stage (Micos) _ 0.05 µm nominal precision


CMOS camera 1280×1024 pixel 256 gray levels 376 frames/sec (Mikrotron MC1310)

Emulsion Plate

XY stage (Micos) 0.1 µm nominal precision



Illumination system, objective (Oil 50× NA 0.85) and optical tube (Nikon)





16 Tomographic Images taken through 44-micron Emulsion Layer

Movable stage

What the microscope CCD sees in one film..


Data (images) processing and motion control flow in the ESS





Kinematical measurements

1. Momentum measurement by multiple coulomb scattering



0.8 GeV/c pion : P = 0.79 (GeV/c), dP/P = 11%

1/p& 0.8GeV/c


1.5GeV/c pion : P = 1.53(GeV/c), dP/P = 16%

1/P(GeV/c)-1

Exposure of an Emulsion Cloud Chamber (ECC) to 400 MeV/u Carbon ions


ECC structure: 219 emulsion films (300 µm thick) and 219 Lexan plates 1 mm thick (73 consecutive "cells")


Lexan: $\rho = 1.15 \text{ g/cm}^3$ and electron density = $3.6 \times 10^{23}/\text{cm}^3$, e.g. Water $3.3 \times 10^{23}/\text{cm}^3$




Carbon exposure at HIMAC (NIRS-Chiba)












Vertex reconstruction



Charge separation






o Co/2 De Lellis et al., Journal of Instrumentation 2 (2007) P06004 13

Charge changing cross-section

130

G. De Lellis et al. / Nuclear Physics A 853 (2011) 124-134



G. De Lellis et al. / Nuclear Physics A 853 (2011) 124–134





Fig. 4. Distribution of the sum of the absolute values of the electrical charge of the particles associated to multi-provertices.

Fig. 2. Fraction of the remnant Carbon beam as a function of the traversed ECC material.


$$\sigma_{tot} = (18420 \pm 380_{stat} \pm 1840_{sys}) \text{ mbarn}$$

$$\sigma(\Delta z = 1) = (2510 \pm 140_{stat} \pm 250_{sys})$$
 mbarn $\sigma(\Delta z = 2) = (1170 \pm 90_{stat} \pm 120_{sys})$ mbarn $\sigma(\Delta z = 3) = (1460 \pm 105_{stat} \pm 150_{sys})$ mbarn $\sigma(\Delta z = 4) = (7510 \pm 240_{stat} \pm 750_{sys})$ mbarn


Emulsions in FIRST: Schematic view




Bricks



Brick structure: 61 films and 55 lead plates




Track slopes expected in emulsion




 $tg(\alpha) > 1 \rightarrow \alpha > 45$ degrees Typical values in high energy physics: $tg(\alpha) < 0.6$

Dedicated study on Large angle track detection




Brick 1: preliminary results



Scanning of 4 emulsion films (slope ≤ 2)


Tracks angular distribution of films 61÷58

Peak signal selection




Scanning of brick 2

• 3 most upstream films




brick rotated by 180 degrees w.r.t. the other one


Brick 2: preliminary results




Scanning of 3 emulsion films (slope ≤ 2)

Tracks angular distribution of films 61÷59





Cosmic-ray angular distribution



Conclusions

- Scanning of very large angle tracks is now routinely done
- First look at the films exposed at GSI in Summer 2011
- To be done:
- Define tracks by the first 6 consecutive films
- Analysis of the films interleaved with lead to measure the properties of these tracks and try the particle identification (proton/alpha)