Precision Timing Calorimetry for High Energy Physics



Adi Bornheim

Dustin Anderson, Cristian Pena, Artur Apresyan, Maria Spiropulu, Javier Duarte, Si Xie, Anatoly Ronzhin

Caltech, FNAL



Precision Timing Calorimetry

4D event reconstruction opens new territory in difficult experimental environments.

78 pp collision LHC bunch crossing

0.11ns 0.14ns - 0.02ns


11ns

0.8ns



Challenges at HL-LHC


- Large samples needed to fully exploit LHC, goal : collect x10 more
 - <PU> ≈ 140 at HL-LHC → 50nb/sec , collect 3000 fb⁻¹
- Some key signatures at HL-LHC
 - Higgs VBF and $W_L W_L$ scattering with forward jets, vertex identification for $H \rightarrow \gamma \gamma$
 - Searches in final states with MET from LSP
 - Precision studies of new physics which may be discovered at LHC








Precision timing at HL-LHC


- Target resolution of O (20-30 psec)
- Allows reconstruction of H→γγ vertex and ~x10 pileup suppression
- Applications of timing information:
- <u>Object level</u> : (e.g. identify forward PU jets for VBF Higgs, WW scattering)
- <u>Hit level</u> : (e.g. timing-based cluster cleaning)
- <u>Event level</u> (hard scatter vertex reconstruction, e.g. for $H \rightarrow \gamma \gamma$)
- <u>Separate</u> spatially overlapping vertices that originate at different times












Photon Traces in LYSO Crystal


- For high energy showers in high light yield crystals, number of scintillation light yield is very large (>10⁵ / GeV).
- Photon detection at one location in the crystal will be an averaged transit time spectrum





Scintillation Light Time Spectrum


- Scintillating crystals get often classified in fast and slow by their light output decay constants. This is often 10s of ns – PWO, LYSO : ~40 ns.
- Timing information is extracted from the leading edge of the signal the rise time of the light output is important.
- > LYSO scintillation light properties :
 - > Light output rise time $t_R < 75$ ps, 35000 photons/MeV, $t_D = 33$ ns.
 - See : S Seifert, J H L Steenbergen, H T van Dam and D R Schaart, 2012 JINST 7 P09004. doi:10.1088/1748-0221/7/09/P09004






Optical Transit Time Spread

- Effect of the scintillation photon arrival at the photo detector we refer to as Optical Transit Time Spread.
- Experimental program to explore ultimate timing resolution, in particular the impact of the optical transit time spread.



Time evolution of a shower from photon in CMS ECAL PbWO crystal (25 cm long).






Shower Fluctuations



- Utilizing the precise time of arrival requires association with a precise spatial information.
- Shower depth fluctuation may get partly compensated by light propagation to the photo detector.
- Shower fluctuations in particular conversion depth for photons need to be considered, relevant dimensions R_M and X₀ are order cm (30 ps).
- Shown here : GEANT hit time spectrum in a solid PbWO crystal matrix, 100 GeV photons.





շ(է₁ - է₂) [ns]

10-1

Timing Performance of CMS ECAL

Large PbWO crystal calorimeter. **Results from pp collision data at LHC :**

> Electron showers from $Z \rightarrow ee$ decay Δt_{TOF} : ~270 ps, single channel : ~190 ps, without path length correction : ~380 ps

 $\sigma(t_1-t_2)[ns]$

10

10²

 10^{3}

 A_{eff} / σ_n

- Constant term of resolution : ~20 ps in test beam, ~70 ps in situ (same clock).
- Studies on jet timing vertex resolution suggest very promising performance.

10

10

 10^{2}

 $\sigma(t_1 - t_2) = \frac{N}{\Lambda - t_2} \oplus \sqrt{2} \overline{C}$

 $\overline{C} = 0.020 \pm 0.004 \text{ ns}$

ndf = 173 / 169

N = 35.1± 0.2 ns


10

E in EB [GeV]

E in EE [GeV]

CMS 2008


10





Fast timing with Microchannel Plates

- Starting point in exploring precision timing in calorimeters
 - Secondary emitter material as active element in a sandwich type calorimeter
 - First proposed: "On possibility to make a new type of calorimeter: radiation resistant and fast", A. I. Ronzhin et. al, preprint IFVE 90-99, 1990.




- Secondary particles from EM shower are detected by MCP
 - > Signal is proportional to the number of secondaries \rightarrow energy of parent
 - > Most of secondary particles are low energy \rightarrow MCP very efficient
 - > MCP are intrinsically very fast \rightarrow calorimeter with very fast timing




Photek 240 and Photonis MCP-PMT

- TOF time resolution for protons between two MCPs (Photonis vs Photek) found to be ~17 ps.
- MCPs in SEC mode better than 40 ps. 200
- Includes 5 ps from readout (DRS4).
- MCPs serve as our reference timing.



A. Ronzhin et. al. NIM A, Vol 749 p 65-73



25 μm pore size, 60x60mm² sensitive area, rise time~300 ps, SPTR~120 ps,



10 µm pore size, 41mm aperture, PC-MCP distance ~5mm, rise time~60 ps, SPTR~40 ps





From single MIP to full showers with MCPs


See poster presentation :

Fermilab: Sergey Los, <u>Erik Ramberg</u>, Erik Ramberg, CalTech: Artur Apresyan, Si Xie, Maria Spiropulu, U. Of Chicago: Heejong Kim And forthcoming NIM paper (submitted) from the same authors.

- Measurements of shower profile with MCPs as active layer.
- Time resolution as a function of the shower depth : ~13 ps with Photek, <40 ps with Photonis.
- Time resolution among different transverse regions inside a shower : ~30 ps with Photonis
- To appear in NIM.










Experimental setup LYSO timing




Study the effect of scintillation (of LYSO) on time resolution

Minimize the effect of optical transit by using a relatively small LYSO crystal (1.7cm x 1.7cm x 1.7cm cube)




TOF Measurements (1.7 cm³ LYSO)







Increased Light Path Complexity

- R&D for HL-LHC on LYSO/W Shashlik calorimeter.
- Radiation hard in HL-LHC environment up to 3 ab⁻¹,energy resolution of 10%/sqrt(E)⊕1%.
- Resolution performance demonstrated in test beam on a 4x4 matrix.
- Radiation hardness of LYSO tested up to 80% of the required dose.
- Use single Shashlik cell to test timing performance with very complex light path.











Beam Test Setup Shashlik Cell



Adi Bornheim, Elba 2015, Precision Timing Calorimeter




Wavelength Shifter Timing Properties



- Pulse rise time of bare LYSO driven by MCP+DRS4, in Shashlik configuration by timing characteristics of the wavelength shifter.
- Ray tracing simulation of the full optical chain, including LYSO and WLS timing properties reproduces the measurement.
- No additional shaping of the pulse due to the complex light path.



Shashlik Timing Performance




- Performance of solid LYSO cube and LYSO/W scales with the rise time difference due to the WLS.
- Few 10 ps resolution achievable with LYSO based calorimeter, reaching ~32 ps at 32 GeV equivalent signal.

NIM A, Vol 794 (2015) p 7-14



Current Systematic Limits


- Extended measurement up to 150 GeV suggest that the systematic limit is small.
- Fit yields constant term compatible with the reference time resolution of around 15 ps.






Light Sensing Performance

- Differential timing resolution of two MCPs read out with DRS4, illuminated with a 50 ps FWHM laser, DRS4 readout : $\sigma_t = 7.2$ ps.
- Hamamatsu MCP : IRF = 45 ps, TTS = 25 ps.
- Similar performance to MIP in a bare MCP.
- Multiphoton timing of MCP resolution approaches limit of readout chain.








Optimized Pulse Reconstruction

- Extract timing from LYSO signal only from the initial ~2.0 ns of the pulse, fitting a linear function to 10 DRS samples of 200 ps.
- Further significant improvement of the performance down to 42 ps at 32 GeV with a reference time precision of about 20 ps.
- We observe ringing like noise in the MCP signal after a few ns. Not present with SiPMs.





Next Steps

- Crystals and MCPs as sensing elements in a sampling configuration.
- Will look at silicon sensors as well. Sensors suitable for calorimeters may not easily achieve 10 ps resolution – however large number of sensors may compensate this.
- Measure timing performance of 4x4 LYSO/W matrix.
- Move towards a MCP setup which covers entire showers to allow energy measurement.



Summary

- We measured the timing resolution of a LYSO/W Shashlik calorimeter cell to be 42 ps for 32 GeV electrons.
- The timing is extracted from the same scintillation light signal used for the energy measurement.
- We measure shower timing of high energy electrons with commercial MCPs at the level of 15 ps.
- We expect further improvement of the LYSO based measurements with a better understanding of the photo sensor characteristics.
- A large scale calorimeter with a time resolution of a few 10 ps for pile-up mitigation at HL-LHC seems achievable.



Backup



Some Related Poster at this Conference

- "Test beam results of micro channel plates in "ionisation mode" for the detection of single charged particle and electromagnetic showers ", Speaker: Paolo Meridiani
- "Energy and time resolution for a LYSO matrix prototype of the Mu2e experiment", Speaker : Simona Giovannella
- "Development of solar blind UV extended APD for readout of Barium Floride crystals", Speaker : Prof. David Hitlin
- "Test and characterization of SiPMs intended as detector for the MEG high resolution timing counter", Speaker : Marcello Simonetta
- "Fast Timing Detector R&D for the HL-LHC era", Speaker : Dr. Sebastian White
- "State of the art silicon photomultipliers with LSO: Ce codoped Ca scintillators achieve 84ps coincidence time resolution for PET" Speaker : Mr. MYTHRA VARUN NEMALLAPUDI