Super-resolution and looking deep in a volume a melting pot for research groups and small enterprize

Eugenio DelRe

Physics Department, Unversity of Rome "La Sapienza", Italy





eugenio.delre@uniroma1.it - Group & Research Initiatives

Group

Giuseppe Di Domenico, PhD Student Single-photon sources – F005

Fabrizio Di Mei, PhD Student

Scale-free optics and giant response – F008

Davide Pietrangeli, PhD Student
Solitons in gratings - F008
Digital Materials -F008
Rogue Waves - F008

Salvatore Silvestri, PhD Student Super-resolved imaging – IIT

Piergiorgio Caramazza, Laurea Degree Student Mariano Flammini, Laurea Degree Student Mario Ferraro, Laurea Degree Student

Projects

Prin2012 – Rogue Waves
FIRB PHOCOS – Spatial Solitons
ISC – Structured Illumination Imaging



Fundamental limits to far-field (at a distance) imaging using light



The Shackles of Diffraction

Diffraction-limited (details below ~ 200 nm are **lost**)





Fig. 6. A field of 50-nm fluorescent beads, imaged by conventional microscopy (a), conventional microscopy plus filtering (b), linear structured illumination (c), and saturated structured illumination using illumination pulses with 5.3 mJ/cm² energy density, taking into account three harmonic orders in the processing (d). Because no scanning is necessary, a wide field can be imaged simultaneously.


Super-resolution to below 50 nm (through nonlinear microscopy)

Imaging using light


Propagation/Optics




Detector Emitter



Imaging using light




Imaging using light




B. Huang, M. Bates, and X. Zhuang, Annu. Rev. Biochem. 2009. 78:993–1016 (2009)






Our "vision" of the world




Science slips beyond conventional optics...



Seeings things below the diffraction/aperture limit at a distance




Combining imaging with prior knowledge





Combining imaging with prior knowledge



Combining imaging with «prior» knowledge – Super-resolution!

But how?

This is a *multiphoton*, i.e., **NONLINEAR**, process!

Not a single image, but a <u>cumulative</u> effect of many images.

One photon $-\lambda$

Two photon – $\lambda/2$

....











Imaging through fluorescence




B. Huang, M. Bates, and X. Zhuang, Annu. Rev. Biochem. 2009. 78:993–1016 (2009)

Super-resolution imaging




Stimulated emission depletion microscopy



Snap-shots of synaptic vescicles

Super-resolution imaging




B. Huang, M. Bates, and X. Zhuang, Annu. Rev. Biochem. 2009. 78:993–1016 (2009)


Super-resolution imaging

Excitation

SIM – Structured Illumination microscopy




Fluorescent




B. Huang, M. Bates, and X. Zhuang, Annu. Rev. Biochem. 2009. 78:993–1016 (2009)




Increasing resolution by *superimposing known* images on *unknown* objects





Digital Age revival of an age-long secret...



The moiré pattern








We DO NOT see the two lines

We DO "see" the two lines

Is resolution really increased?



TIRF image of a living S2 cell expressing GFP-tubulin.



2D SIM

R. Heintzmann and M. G. L. Gustafsson, Nature Photonics 3, 362(2009)



Increasing resolution by *superimposing known* images on *unknown* objects


Structured Illumination Microscopy








Patterns






Setup





Sample results







Particular of BPAE cells



Sample results





Particular of Human Skin






Scale-free optics




LETTERS

Scale-free optics and diffractionless waves in nanodisordered ferroelectrics

E. DelRe^{1*}, E. Spinozzi¹, A. J. Agranat² and C. Conti³




Comparing standard optics and scale-free optics





$$-\frac{\nabla^2 E}{E} + \left(\frac{L}{\lambda}\right)^2 \left(\frac{\nabla |E|^2}{2|E|^2}\right)^2 = k^2$$



Observation of diffraction cancellation for nonparaxial beams in the scale-free-optics regime

F. Di Mei, ^{1,2} D. Pierangeli, ¹ J. Parravicini, ¹ C. Conti, ^{3,1} A. J. Agranat, ⁴ and E. DelRe^{1,3,*}


Challenge

Propagate a non-paraxial subwavelength-sized beam in a volume along macroscopic distances for imaging




Diffraction

 $\mbox{d} >> \lambda \label{eq:lambda}$ Small angular aperture




 $\mbox{d} \approx \lambda \label{eq:large_large}$ Large angular aperture







?iffraction

? angular aperture




? angular aperture







Material

Dipolar glass forming inside the KTN:Li crystal





Compositional disorder





deep-in-band Cu impurities→Huge non-linear EO response


E. DelRe et al. Nat. Photon. 5, 39 (2011)

A.A. Bokov, Z. Ye Journal of Material Science 41 31-52 (2006)




Subwavelength anti-diffracting beams propagating over more than 1,000 Rayleigh lengths

Eugenio DelRe^{1*}, Fabrizio Di Mei^{1,2}, Jacopo Parravicini¹, Gianbattista Parravicini³, Aharon J. Agranat⁴ and Claudio Conti^{1,5}









50 μm

Massive beam propagation (L/ λ =1.1)





Subwavelength massive beam propagation (L/ λ =1.1)






Structured illumination 3D scan



Idea: structured light with periodic on-off pattern


Experimental setup





Vertical strip pattern


Strips distance versus surface displacement




Surface displacement produces relative distance modifications between pattern strips coming from two different directions




Test surface and acquisitions




Elastic membrane deformed with spherical profile in range 0 to 3 mm







Processing algorithm




Peaks detection and identification (low and high)

Peak alignment respect frame n

stripes distances processing

Displacement processing



Processing is performed using parallel algorithms to improve efficiency

