

Recent measurements using monoenergetic and thermal neutrons at the National Physical Laboratory

Nigel Hawkes, Catalin Matei, Neil Roberts, Graeme Taylor, David Thomas

Bushy House





NPL was founded in 1900 in a former royal residence

New laboratories





Completed 2010

Neutron facility






Chadwick Building


3.5 MV Van de Graaff accelerator





Showing ion source, pulser and accelerator tube


Experimental area



Showing low scatter area for monoenergetic & radioisotope neutrons (left side) and thermal pile (right)

Low-scatter area





At least 6 m from the neutron source to the shield walls & real floor. Low-mass detector supports & walkways.

Neutron fluence

National Physical Laboratory

- Measured using a Long Counter (BF₃ counter inside a cylindrical moderator)
- The efficiency, about 11 counts per (neutron per cm²), varies relatively little with energy
- Efficiency established by (e.g.) using neutron sources with accurately known output (measured in NPL Mn bath)



Monoenergetic neutrons



Typically used for:

- Calibrating external neutron standards
- Testing and characterising novel detectors
- Nuclear data

Specific example:

- •Nuclear modelling relies on accurate and precise data.
- ■Calculations for Generation IV nuclear reactors require
 ²⁴²Pu fission cross section to 3 5 percent.
- ■But currently the uncertainty is 19 21 percent.
- •This is currently being addressed in an IRMM / NPL / PTB collaboration as part of the EMRP Metrofission project


²⁴²Pu target in ionization chamber NPL






Experimental setup





²⁴²Pu (n,f) cross section






Thermal Pile







Graphite block about 2.8 m long by 1.4 m wide by 1.6 m high

Thermal Pile



Central cavity:

- Isotropic field
- Up to 2×10⁷ cm⁻² s⁻¹
- Diameter ~12 cm
- Cd ratio ~33


Thermal column:


Beam geometry

Up to 4×10⁴ cm⁻² s⁻¹

Diameter ~30 cm Boronated


Cd ratio ~6.5






Central cavity – fluence rate (at max setting) & Cd ratio






MCNP simulation

- Behaviour of Thermal Pile modelled using MCNP
- Neutrons that enter the central cavity commonly do so multiple times (2.3 on average)







Measuring the fluence delivered



Gold foil activation followed by off-line β –counting in a 4π low background β –counter (or β – γ counter).





Conclusions



- The NPL Neutron Metrology Group has a 3.5 MV Van de Graaff accelerator and a range of experimental facilities.
- Well-characterised neutron fields can be produced in the fast and thermal energy regions.
- These can be used for calibrations, activations, reactor instrument testing, and cross section measurements.